# ESK+ Stanceflex

## Instructions for Use

ESKSFP30  
ESKSFP30PYR

SMARTIPSF30  
SMARTIPSF30PYR

ESKSF30  
ESKSFP30PYR

EN	Instructions for Use	2
DA	Brugsanvisning	20
NO	Bruksanvisning	38
FI	Käyttöohjeet	56
SV	Bruksanvisning	74
RU	Инструкция протезиста	92
ZH	使用说明	110
AR	إرشادات الاستخدام	128

Contents .....	2
1 Description and Intended Purpose .....	3
2 Safety Information .....	5
3 Construction .....	6
4 Function.....	7
5 Maintenance .....	7
6 Limitations on Use .....	8
7 Bench Alignment.....	9
8 Alignment Advice and Stability Adjustment.....	10
9 Fitting Advice.....	11
9.1 Changing Stanceflex Rubbers .....	11
9.2 Front Stop Adjustment - Adjustable Trunnion .....	12
9.3 Further Adjustments.....	13
10 Assembly Instructions .....	14
10.1 Changing Front Stop Buffer.....	14
10.2 Changing Shin Stops.....	14
10.3 Removing ESK+ Mechanism .....	15
10.4 Replacement of ESK+ Mechanism .....	16
11 Technical Data .....	17
12 Ordering Information .....	18

---

# 1 Description and Intended Purpose

These Instructions for Use are intended for use by the practitioner and user unless otherwise stated.

The term device is used throughout these instructions for use to refer to Stanceflex ESK+ products.

Please read and ensure you understand all instructions for use, in particular all safety information and maintenance instructions.

## Application

This device is to be used exclusively as part of a lower limb prosthesis.

This device is for a single user.

The device is a uniaxial weight activated stabilized knee with integrated stance flexion element and is intended to improve safety and reduce compensation.

The shin can be fitted with a PSPC (Pneumatic Swing Phase Control) or SmartIP (Smart Intelligent Prosthesis) swing phase control or supplied without a cylinder.

A rubber element at the rear of the knee chassis provides extra stance flexion.

## Features

	Manually adjusted PSPC	Microprocessor controlled SmartIP
Stance flexion during weight bearing	✓	✓
Tough carbon fibre frame	✓	✓
140 degree knee flexion	✓	✓
Adjustable sensitivity of knee stability	✓	✓
Independent manually adjustable swing controls for flexion and extension	✓	
Microprocessor controlled swing control in flexion with manually adjustable terminal impact damping in extension		✓

## Activity Level

This device is for Activity Level 2, 3 and 4 users (weight limits apply, see *Technical Data*).

This device might not be suitable for Activity Level 1 users or for use in competitive sports events. These types of users might be better served by a specially designed prosthesis that is optimized for their needs.

### Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

### Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low-level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

### Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

### Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

## Clinical Benefits

- Weight-activated stance control and adjustable swing control provide smooth, safe gait with a free knee.
- Stanceflex improves shock absorption.
- Microprocessor swing control allows optimized gait characteristics across a range of walking speeds. (SmartIP only.)

## Contraindications

This device is not suitable for:

- Users with extreme muscular weakness, contractures that cannot be correctly accommodated, or proprioceptive dysfunction including poor balance.
- Users with contralateral joint instabilities or pathology.
- Users with complicated conditions involving multiple disabilities.
- Users, including bilateral transfemoral amputees, who require a yielding function from the knee to descend stairs or to sit. This can result in damage of the brake mechanism that can adversely affect stability and durability.

**Users must be given gait training before using this device.**

**Make sure that the user has understood all instructions for use, drawing particular attention to all maintenance and safety information sections.**

## 2 Safety Information

**This warning symbol highlights important safety information.**

**Beware of finger trap hazard all times.**

**Any changes in the performance or function of the limb e.g. instability, double-action, restricted movement, non-smooth motion or unusual noises should be immediately reported to your service provider.**

**Always use a hand rail when descending stairs and at any other time if available.**

**Any excessive changes in the heel height after finalization of alignment may adversely affect limb function.**

**Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.**

**Care should be taken when carrying heavy loads as this may adversely affect the stability of the device.**

**The user should be advised to contact their practitioner if their condition changes.**

**Avoid exposure to extreme heat and/or cold.**

**Assembly, maintenance and repair of the device must only be carried out by a suitably qualified clinician.**


**The device is not intended for use when immersed in water or as a shower prosthesis. If the limb comes into contact with water wipe dry immediately.**

# 3 Construction

## Principal Parts

- Carrier Assembly Carbon fibre composite & stainless steel
- Chassis Assembly Aluminium alloy, stainless steel, engineering thermoplastic, fibre reinforced (asbestos free) friction material.
- Distal Dome Aluminium alloy
- Distal Pyramid Titanium with stainless steel bolt
- Distal Tube Clamp Aluminium alloy
- Patella Engineering thermoplastic
- Stanceflex Rubber Natural rubber
- Stanceflex strap Kevlar
- PSCP Aluminium alloy, stainless steel, titanium, engineering thermoplastic
- SmartIP Control Cylinder Aluminium alloy, stainless steel, Titanium, engineering thermoplastic  
Electronic components, batteries

---

## 4 Function

The knee is a uniaxial design comprising a knee pivot connecting a chassis to a carrier. A weight activated stabilizing knee mechanism is housed within the chassis. The knee chassis has a pivotally connected top plate which is resiliently supported by a rubber block, this provides a stance flexion action under load and absorbs shock loads at heel strike.

A knee control cylinder sits within this frame assembly, fitted between the chassis posterior pivot and the carrier distal pivot.

The PSPC cylinder has manually adjustable valves which give independent control of flexion and extension resistance during swing phase.

The SmartIP control cylinder uses sensors inside the cylinder to monitor walking speed and automatically adjusts flexion resistance accordingly. Extension control is via fixed adjustment of a terminal impact damper.

---

## 5 Maintenance

Visually check the device regularly.

Report any changes in performance of this device to the practitioner/service provider e.g. unusual noises, increased stiffness or restricted/excessive rotation, significant wear or excessive discoloration from long term exposure to UV.

Inform the practitioner/service provider of any changes in body weight and/or activity level.

### Cleaning

Use a damp cloth and mild soap to clean the outside surfaces. DO NOT use aggressive cleansers.

***The remaining instructions in this section are for practitioner use only.***

This maintenance must be carried out only by competent personnel (practitioner or suitable trained technician).

To maintain optimum performance, a periodic inspection for unacceptable noise, play and stability is recommended. The period between inspections is determined by factors such as the activity of the user, frequency of usage and the stability setting of the knee.

As a guide, if the setting is 6 clicks or above and the user averages 600,000 steps per annum, an annual inspection is recommended.

If the stability setting is less than 6 clicks and/or the user averages fewer than the average number of steps, the period in between inspections may be longer.

The use of a Limb Activity Monitor is recommended to determine step count.

We recommend adjusting the knee to the individuals' stability requirements to ensure optimal stability and durability. See adjustment guide/flow diagram.

At these inspections:

- Check for visual defects that may affect proper function.
- Check for M-L play of chassis and excess compression of the Stanceflex rubber. Replace bearings and Stanceflex rubber if necessary.
- Check front stop and replace if necessary.
- Check shin stops and replace if necessary.
- Dis-assemble, clean and reassemble knee mechanism.

Advise the user that a regular visual check of the device is recommended and any changes in performance of this device must be reported to the practitioner. Changes in performance may include: increase or decrease in knee stiffness, instability, reduced knee support (free movement), any unusual noise.

If this device is used for extreme activity, the maintenance level and interval should be reviewed and if required advice and technical support sought to plan a new maintenance schedule dependent upon the frequency and nature of the activity. This should be determined by a local risk assessment carried out by a suitably qualified individual.

### **Storage and Handling**

When storing for prolonged periods, the device should be kept dry and within a temperature range of -15 °C to 50 °C. (5 °F to 122 °F)

Use product packaging supplied.

---

## 6 Limitations on Use

### **Intended life:**

A local risk assessment should be carried out based upon activity and usage.

### **Lifting Loads**

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

### **Environment:**

Avoid exposing the device to corrosive elements such as water, acids and other liquids. Also avoid abrasive environments such as those containing sand for example as these may promote premature wear.

Exclusively for use between -15 °C and 50 °C  
(5 °F to 122 °F)

Suitable for outdoor use


# 7 Bench Alignment

*The instructions in this section are for practitioner use only.*

## Walking Trials

### Spring Stack Adjuster

Fully extend limb to access the adjuster.


Rotate the **SPRING STACK ADJUSTER** to achieve the required stability for the individual amputee.

Note... Audible clicks should be heard upon adjustment, and the adjuster is self locking.

**Tighten the SPRING STACK to reduce stability.**


**Release the SPRING STACK to increase stability**  
- See above.

If required, the adjuster may be reset to the factory setting. This can be achieved by reducing stability by six audible clicks from the maximum stability setting.

# 8 Alignment Advice and Stability Adjustment


The instructions in this section are for practitioner use only.

# 9 Fitting Advice



*The instructions in this section are for practitioner use only.*

## 9.1 Changing Stanceflex Rubbers

Static





Knee Flexion Following Heel Strike

## 9.2 Front Stop Adjustment - Adjustable Trunnion

*The instructions in this section are for practitioner use only.*

### Fitted to the PSPC

### Front Stop Adjustment - Non-Adjustable Intelligent Prosthesis

(As fitted to older SmartIP models.)

No adjustment is necessary when using these cylinders. They require the front stop to be fully tightened as shown in Stage 1.

## 9.3 Further Adjustments

*The instructions in this section are for practitioner use only.*

Symptom	Solution
Knee does not stabilise	Increase sensitivity
	Shift socket Anteriorly
	Replace stabilizing mechanism if trailing edge of brake is damaged
	Check Stanceflex rubber and replace if necessary (See minimum thickness in compression in <i>Knee Flexion Following Heel Strike</i> section)
Knee does not release	Decrease sensitivity
	Check correct adjustment of piston rod length
	Shift socket posteriorly
Heel rise too high during swing	See control cylinder IFU
Heel rise too low during swing	See control cylinder IFU
User experiences terminal impact	See control cylinder IFU
Device does not fully extend during swing	See control cylinder IFU



If after following these instructions you still have problems with the function please contact the sales team in your area for advice.

# 10 Assembly Instructions


*The instructions in this section are for practitioner use only.*

## 10.1 Changing Front Stop Buffer


When fitting a new front stop, always re-adjust, and if applicable check length of cylinder if fitted with an adjustable trunnion. (See previous page.)

## 10.2 Changing Shin Stops









**1** Remove shin stops and remainder of adhesive. Clean and lightly abrade mating surfaces.

**2** Apply Loctite 424 on mating surfaces. Press shin stops into place.









## 10.3 Removing ESK+ Mechanism

*The instructions in this section are for practitioner use only.*

<p><b>1</b> Remove all preload and record number of clicks</p>  <p>Underside view of chassis</p>	<p><b>2</b></p> 	<p><b>3</b></p> 
<p><b>4</b></p> 	<p><b>5</b></p> 	<p><b>6</b></p> 
<p><b>7</b></p>  <p>Oscillate during removal</p> <p>Brake band material may be wiped clean with paper, do not use solvents.</p>		<p><b>8</b></p>  <p>Clean bore with suitable solvent</p>

## 10.4 Replacement of ESK+ Mechanism

*The instructions in this section are for practitioner use only.*


<p><b>1</b></p>  <p>Offer up side plate to align slot in spindle (2) with dowel (1) but do not fit at this stage.</p>	<p><b>2</b></p>  <p>Screw spindle to chassis.</p>	
<p><b>3</b></p> <p>Oscillate during assembly</p>  <p><b>ORIENTATION IS IMPORTANT!</b> Assemble flat sides and lettering (RH, LH) facing inwards.</p>	<p><b>4</b></p>  <p>LOCTITE 222</p>	<p><b>5</b></p>  <p>STP Oil Treatment</p>
<p><b>6</b></p> 	<p><b>7</b></p> 	<p><b>8</b></p> <p>Re-apply preload with recorded number of clicks.</p> <p>Underside view of chassis.</p>  <p>2.5</p>


# 11 Technical Data

Operating and Storage Temperature Range:	-15 °C to 50 °C (5 °F to 122 °F)
Component Weight:	
SmartIP	1.27 kg (2 lb 12.8 ozs)
PSPC	1.13 kg (2 lb 8 ozs)
Activity Level:	2-3 (4)
Maximum User Weight:	
Activity Levels 2-3	125 kg (275 lb)
Activity Level 4	100 kg (220 lb 7.4 ozs)
Maximum Flexion Angle:	140°
Proximal Alignment Attachment:	Single Bolt (Blatchford)
Distal Alignment Attachment:	30 mm Tube Clamp or Male Pyramid (Blatchford)
Build Height:	See Diagrams Below.

## Fitting Length

**ESK+ Stanceflex  
with Distal Pyramid**

**ESK+ Stanceflex  
with 30 mm Distal Tube Clamp**

# 12 Ordering Information

## Order Example

	ESK+ with Stanceflex							
	30 mm Distal Tube Clamp				Distal Pyramid			
Pneumatic Cylinder	ESK	SF	P	30	ESK	SF	P	PYR
	ESK+	Stanceflex	PSPC	Distal Tube Clamp 30 mm	ESK+	Stanceflex	PSPC	Distal Pyramid
SmartIP Knee	SMARTIP	SF		30	SMARTIP	SF		PYR
	SmartIP Knee	Stanceflex		Distal Tube Clamp 30 mm	SmartIP Knee	Stanceflex		Distal Pyramid
No Control Cylinder Fitted	ESK	SF		30	ESK	SF		PYR
	ESK+	Stanceflex		Distal Tube Clamp 30 mm	ESK+	Stanceflex		Distal Pyramid

## Replacement Parts

Item	Part. No.
Stanceflex Knee Assembly	239149
ESK Mechanism Alignment Pin	234914
Front Stop Pad	234903
ESK+ Stabilizing Mechanism Kit	239949
Shin Stop	233810

## Stanceflex Rubber

Item	Part. No.
25° Pink (Soft)	231608
35° White	231606
45° Orange (Fitted as Standard)	231607
70° Blue (Hard)	231609

## Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

## CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: [www.blatchford.co.uk](http://www.blatchford.co.uk)

Medical Device

Single Patient – multiple use

## Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

## Warranty

This device is warranted for 24 months, the brake band for 12 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

See Blatchford website for the current full warranty statement.

## Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

## Environmental Aspects

This product contains hydraulic oil, mixed metals and plastics. Where possible the components should be recycled in accordance with local waste handling regulations.

## Retaining the Packaging Label

The practitioner is advised to keep the packaging label as a record of the device supplied.

## Manufacturer's Registered Address

Blatchford Products Limited, Lister Road, Basingstoke RG22 4AH, UK.

## Trademark Acknowledgements

Blatchford is a registered trademark of Blatchford Products Limited.

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**Blatchford Europe GmbH**
**Endolite India Ltd.**