Silcare Active Cushion Liner

Silcare Active Locking Liner

Instructions for Use

SATTCP22-SATTCP40
SATTLP22-SATTLP40

EN	Instructions for Use	2
DA	Brugsanvisning	13
NO	Bruksanvisning	24
FI	Käyttöohjeet	35
SV	Bruksanvisning	46
RU	Инструкция протезиста	57
ZH	使用说明	68

Contents

Contents	2
1 Description and Intended Purpose	3
2 Safety Information	5
3 Construction	6
4 Function.....	6
5 Maintenance	6
5.1 Cleaning the Device	7
5.2 Cleaning the Residual Limb.....	7
6 Limitations on Use	8
7 Choosing the Correct Size.....	8
8 Trimming the Device.....	8
9 Donning the Device	9
10 Fitting Advice.....	9
11 Technical Data	10
12 Ordering Information	11

1 Description and Intended Purpose

These instructions are for use by the practitioner.

Please read these instructions carefully before fitting the device.

The term *device* is used throughout these instructions for use to refer to Silcare Active Cushion/Locking Liner.

Ensure that the user understands all instructions for use, drawing particular attention to the maintenance and safety information.

Application

This device is an interface component for use only as part of a lower limb prosthesis.

Intended for single user.

It is intended for moderate to high activity users with well covered residual limbs.

The Blatchford Silcare Active silicone liner provides a cushioned socket interface and is manufactured from biocompatible materials. Please read these instructions carefully before fitting the device.

To optimize comfort ensure that the user is instructed in the correct way to handle and don/doff the device. See Section 9 *Donning the Device*. Also ensure that they are made aware of how to maintain the device and to keep it clean for maintaining hygiene as set out in these instructions. See Section 5 *Maintenance*.

Activity Level

The device is recommended for Activity Level 3 and 4 users, although may be suitable for all four activity levels.

The silicone used in the Silcare Walk liners is softer and more flexible than that used in Silcare Active liners. The softer silicone offers more comfort and flexibility in all directions but may feel less connected or secure at higher activity levels. Silcare Walk liners are therefore generally recommended for the less active users (Activity Level 1 and 2).

The silicone used in Silcare Active liners is firmer and slightly less flexible, offering a more secure interface at higher activity levels. Silcare Active liners are therefore generally recommended for the more active users (Activity Level 3 and 4).

Note that some high activity users with poor residual limb tissue coverage and/or high sensitivity may prefer the comfort of the softer silicone Silcare Walk liner to the firmer Silcare Active liner. Conversely, some low activity users with good residual limb tissue coverage and/or low sensitivity may prefer the security of the firmer Silcare Active liner to the softer Silcare Walk liner.

Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low-level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

Contraindications

- Conical residual limbs - this device may not be optimally contoured to fit conical residual limbs
- Users with poor hand or cognitive function might find donning and cleaning difficult.
- Poor hygiene
- Users whose suspension sleeve causes a bulk around their knee may be better suited with alternative type liner

Clinical Benefits

- Provides cushioning for the residual limb in the socket
- Distributes in-socket pressure more evenly, compared other materials and alternative cushioning solutions
- The locking liner provides a means of suspension

2 Safety Information

This warning symbol highlights important safety information.

The user should be advised to contact their practitioner if their condition changes.

Any deterioration in residuum condition or any change in sensation should be reported to the practitioner. If the residual limb has tissue damage, the use of the device should be stopped.

Make sure that any damaged skin or open wound is properly and suitably dressed to prevent direct contact with the device.

Wearers with sensitive skin, diabetics and vascular cases should be extra vigilant and may need to apply lubricant to sensitive areas. We recommend a routine visual check and if required the user should consult their medical practitioner.

For other medical conditions the user should follow the advice and recommendation of a physician or a medical practitioner regarding skin care.

Do not use alcohol sprays, household cleaners, or abrasives. These cleaning materials could damage the device and irritate skin.

Do not pull or stretch the fabric. Fingernails, sharp jewelry, and the locking pin can tear the fabric. If the fabric is torn, stop using it and contact a Blatchford sales representative.

Sockets with sharp proximal edges can tear the device.

Take care when handling the device to avoid possible contamination with materials such as fibreglass which will stick to the device and cause skin irritation.

When donning a sock, clothing and the prosthetic limb, be aware that the device can build up a static charge.

To avoid danger of suffocation, keep the device away from babies and children.

Keep the device away from direct heat sources.

Only use the device in combination with corrosion resistant components.


Do not overtighten the locking pin.

3 Construction

Principal Parts

- Fabric (polyamide and lycra)
- Main Body (silicone)
- Distal Cap (silicone)
- Umbrella (nylon)

4 Function

The device, used in connection with an air-tight suspension sleeve (cushion liner only) and corrosion resistant components, provides control, secure connection and cushioning between the limb and the socket.

5 Maintenance

Advise users to report the following to their practitioner:

- Tears in the fabric or the silicone
- Tissue damage on the residual limb
- Changes in either body weight or activity level
- Deterioration/changes to the residual limb
- Changes in the performance of the device

Note... Perspiration discolors some socket materials.

5.1 Cleaning the Device

Wash the inside of the device daily to avoid any build up of bacteria.

 **Take care handling the device when it is inside out to avoid picking up dust, grit and other contamination which could cause irritation to the skin.**

 **Dry thoroughly before use.**

 **Do not tumble dry.**

Washing by Hand

1. Invert the device so that its silicone side points externally.
2. Clean the silicone and the distal cap with a solution of water and unperfumed, pH balanced soap.
3. Invert the device so that its silicone side points internally.
4. Clean the valve and carefully clean the distal cap (only locking type).
5. Rinse the device with clean water to remove all residues.
6. Either pat dry the device with a lint-free cloth, or leave it to air dry. Take care when drying and handling the device.

Note... Always dry the device with the silicone side of the device pointing internally. Otherwise, the fabric can stretch and become distorted.

Machine Wash

Suitable for machine washing at 30 °C.

 **Do not invert the device for machine washing.**

5.2 Cleaning the Residual Limb

Inspect residual limb before and after prosthetic limb use or at least daily.

 **Any deterioration in residuum condition should be reported to the practitioner. If the residual limb has tissue damage, the use of the device should be stopped.**

1. Clean the skin daily with an unperfumed, pH balanced soap.
2. Rinse the skin with clean water to remove all residues.
3. Apply lotion to dry skin, as recommended by a practitioner.

 **Make sure that any damaged skin or open wound is properly and suitably dressed to prevent direct contact with the device.**

6 Limitations on Use

Intended Life

A local risk assessment should be carried out based upon activity and usage.

Environment

Avoid exposing the device to corrosive elements such as acids, industrial detergents, bleach or chlorine. The use of creams or lotions with this device should be used with caution as these can cause the device to soften and overstretch or distort.

Keep away from sharp objects (such as jewellery, fingernails).

Exclusively for use between -15 °C and 50 °C (5 °F to 122 °F).


Suitable for shower use

7 Choosing the Correct Size

It is important that the correct size is specified to ensure proper fit and comfort in use.

1. Measure the circumference of the residual limb 4 cm up from its distal end.
2. If the circumference of the residual limb matches one of the sizes in the table, choose the next size down.
3. If the circumference of the residual limb is between one of the sizes in the table, choose the smaller size.

Circumference Measured (cm)	22.2-23.5	23.6-25	25.1-26.5	26.6-28	28.1-30	30.1-32	32.1-34	34.1-36	36.1-40	40.1-42
Available Sizes	22	23.5	25	26.5	28	30	32	34	36	40
Valve Size (Locking Liner)	Small			Medium		Large			Extra Large	

8 Trimming the Device

The device may be trimmed to suit as required but never bellow the socket trim line.

We recommend using a rotating disc cutter to achieve a bevelled edge.

Leave enough material on the device so that its proximal edge can cover the edge of the socket.

Trim the proximal edge of the device to suit the shape of the user's socket.

If necessary, trim the proximal edge of the device in a wave pattern to reduce shear stress.

 **Where possible try to avoid abrupt changes of contour and sharp edges which could cause tears in the silicone and fabric.**

9 Donning the Device


 Care must be taken during donning/doffing not to damage the device with fingernails, sharp jewelry or the locking pin.

 Do NOT pull or stretch the device.

1. Roll the device inside out.
2. Make sure the device is positioned correctly to suit the top profile.
3. Roll the device on taking care to expel any trapped air.

After donning the device, allow the user to wear it for 10 minutes. If the user feels numbness, tingling or any unusual sensation within 10 minutes, do the following:

4. Doff the device, and wait until normal sensation returns.
5. Don the device again.
6. If the user feels numbness, tingling or any unusual sensation again, stop using the device.

10 Fitting Advice


Cushion

The device must be worn in conjunction with an air-tight sleeve.


Locking

When using a locking liner the integrated matrix should finish just below fibula head/patella tendon level.

The locking pin must always follow the long axis of the residual limb. (See diagrams.)

The pin should be secured by applying Loctite 243 to the threads of the locking pin. Torque the locking pin to 3 Nm, or tighten it with your fingers by 1/8 to 1/4 of a turn.


11 Technical Data

Principal Materials	polyamide, lycra, silicone, nylon
Shore Hardness	40 shore 00
Component Weight (<i>size 28</i>)	695 g (1 lb 8 oz)
Activity Level	3-4
Size Range (see <i>Choosing the Correct Size</i> section)	22-40 cm
Length (See diagram below)	435 mm
Internal Length (See diagram below)	420 mm
Matrix Length (Locking Type Only)	Approx. 10 cm
Operating and Storage Temperature Range	-15 °C to 50 °C (5 °F to 122 °F)
Distal End Attachment (Locking Type Only)	M10 Requires Locking Pin With a Shoulder*
Locking Pin Shoulder Diameter	13.5-19 mm

*Locking Pin Not Supplied

Dimensions

12 Ordering Information

Order Example

S	A	TT	C	P	26
Silcare	Active	Transtibial	C=Cushion/ L=Locking	Parallel	Size (22-40)*

Available from size 22 to size 40:
SATTCP22 to SATTCP40
SATTLP22 to SATTLP40

* SATTCP23 is for size 23.5. SATTCP26 is for size 26.5.

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: www.blatchford.co.uk

Medical Device

Single Patient – multiple use

Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

Warranty

The device is warranted for 6 months from original date of purchase unless otherwise stated. See Blatchford catalogue for details.

Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

Environmental Aspects

This product is fabricated from silicone rubber and fabric that cannot be easily recycled: please dispose of it responsibly as general waste, as per local waste handling regulations.

Retaining the Packaging Label

You are advised to keep the packaging label as a record of the device supplied.

Trademark Acknowledgments

Silcare and Blatchford are registered trademarks of Blatchford Products Limited.

Manufacturer's Registered Address

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Blatchford Europe GmbH
Endolite India Ltd.