# M36 Adapter Kit

## Instructions for Use

239092

EN	Instructions for Use	2
FR	Instructions d'utilisation	12
DE	Gebrauchsanweisung	22
IT	Istruzioni per l'uso	32
ES	Instrucciones de uso	42
NL	Gebruiksaanwijzing	52
PL	Instrukcja użytkowania	62
PT	Instruções de utilização	72
CS	Návod k použití	82

# Contents

Contents .....	2
1 Description and Intended Purpose .....	3
2 Safety Information .....	4
3 Construction .....	4
4 Function.....	5
5 Maintenance .....	5
6 Limitations on Use .....	6
7 Fitting Advice.....	6
8 Technical Data .....	10
9 Ordering Information .....	10

---

# 1 Description and Intended Purpose

These Instructions for Use are intended for use by the practitioner and user unless otherwise stated.

The term *device* is used throughout these instructions for use to refer to the M36 Adapter Kit.

Please read and ensure you understand all instructions for use, in particular all safety information and maintenance instructions.

## Application

This device connects prosthetic devices and allows for rotational adjustment. The device is to be used exclusively as part of a lower limb prosthesis.

## Features

- When tightened securely the device offers a secure way to connect two parts of the limb.
- The rotational movement offered by the adapter and A-P slide adjustment when fitted to a Blatchford chassis allows for a customized limb adjustment to suit the user.

## Activity Level

This device is suitable for activity levels 1 to 4; weight limits apply, see Section 8 *Technical Data*. Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances and any such decision should be made with sound and thorough justification.

### Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

### Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low level environmental barriers such as, curbs, stairs or uneven surfaces. Typical of the limited community ambulator.

### Activity Level 3

Has the ability or potential for ambulation with variable cadence.

Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

### Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

## Contraindications

There are no known contraindications if used in accordance with these instructions.

## Clinical Benefits

- Provides a connection between the socket and the prosthetic knee.

## 2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Avoid exposure to extreme heat and/or cold.

To reduce the risk of injury due to failure or loosening of the screw connections ensure the bolt threads are cleaned thoroughly before each installation.

The device is designed for prolonged submersion and suitable for immersion in fresh water only. Ensure any use of the device in water complies with the conditions given in Section 6 *Limitations on Use*.

Always apply the specified torque value to the bolts. Never replace the bolt with alternative screws.

Do not use the product if you notice any unusual noises or excessive play in the product.

Ensure the interfacing M36 threaded adaptor is fully engaged onto the thread of the adapter.

The user must not adjust or tamper with the setup of the device.


The user should be advised to contact their practitioner if their condition changes.

## 3 Construction

### Principal Parts

Chassis Bolt	–	Stainless Steel
M36 Proximal Adaptor	–	Aluminum
T-Nut	–	Stainless Steel

### Component Identification

---

## 4 Function

The rotational and slide adjustment allows for a customized limb alignment to suit the user. This device offers a secure way to connect two parts of the limb.

---

## 5 Maintenance

Visually check the device regularly.

Report any changes in performance of this device to the practitioner/service provider e.g.

- Instability
- Excessive play or loss of alignment
- Any unusual noise

Inform the practitioner/service provider of any changes in body weight and/or activity level.

### Cleaning

Use a damp cloth and mild soap to clean outside surfaces, do not use aggressive cleansers.

***The remaining instructions in this section are for practitioner use only.***

This maintenance must be carried out only by competent personnel (practitioner or suitable trained technician).

The following routine maintenance is to be carried out at least annually:

- Check the chassis bolt for tightness. If loose, remove and clean, then apply Loctite and tighten to the correct torque setting.
- Check for defects that could affect proper function.
- Check any mating chassis component screws are tight (see diagrams in Section 7 *Fitting Advice*).

Ensure the user has read and understood all safety and user-level maintenance information.

Advise the user that a regular visual check of the device is recommended and any changes in performance and signs of wear that may affect function should be reported to their service provider.

Advise the user to inform the practitioner/service provider of any changes in body weight and/or activity level.

If this device is used for extreme activity, the maintenance level and interval should be reviewed and if required advice and technical support sought to plan a new maintenance schedule dependent upon the frequency and nature of the activity. This should be determined by a local risk assessment carried out by a suitably qualified individual.

### Storage and Handling

When storing for prolonged periods ensure the product is free from moisture and stored at room temperature.

---

## 6 Limitations on Use

### Intended Life

A local risk assessment should be carried out based upon activity and usage.

### Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

### Environment

The device is suitable for submersion in fresh water only. Light surface corrosion affects neither the function nor the security of this device. However, if it is heavily corroded, stop using and contact your practitioner.

Avoid exposing the device to abrasive or corrosive environments, such as those containing sand, salt or chemically treated water e.g. swimming pools. After any exposure to such environments rinse the device in fresh water and dry thoroughly. Exclusively for use between -15 °C and 50 °C (5 °F and 122 °F).

Suitable for submersion

---

## 7 Fitting Advice

*The following instructions in this section are for practitioner use only.*

Symptom	Solution
A recurring noise occurs between the socket and the knee	The bolt must be tightened to the torque ratings specified in this document. Apply Loctite.
The adapter moves has moved out of position.	Advise the user not to use the device until adjusted, repaired or replaced. Inspect device. If components are worn or damaged, then remove and replace with new components. Reposition as necessary and retighten as per instructions later in this section.

# Linx, Orion3

Linx


938353



Orion3


938374

# Dynamic Pylon




938306

# Mercury




938251

# ESK+ 4 Bolt


938205


# KX06V2

938409

---

## 8 Technical Data

Operating and Storage Temperature Range: -15°C to 50°C (5°F to 122°F)  
Component Weight: 85 g (0 lb 3 oz)


Recommended Activity Level: 1–4

Maximum User Weight: Activity level 4: 100 kg (220 lb)  
Activity Levels 1-3: 125 kg (275 lb)

Range of Adjustment: 360° rotation

---

### Build Height

---

## 9 Ordering Information

Device	
M36 Adapter Kit	239092

---

## Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

## CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: [www.blatchford.co.uk](http://www.blatchford.co.uk)

Medical Device

Single Patient – multiple use

## Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural tests, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

## Warranty

This device is warranted for 24 months. The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses, and exemptions. See the Blatchford website for the current full warranty statement.

## Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

## Environmental Aspects

This product contains mixed metals. It should be recycled where possible in accordance with local waste recycling regulations.

## Retaining the Packaging Label

The practitioner is advised to keep the packaging label as a record of the device supplied.

## Trademark Acknowledgements

Blatchford is a registered trademark of Blatchford Products Limited.

## Manufacturer's Registered Address

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**Blatchford Europe GmbH**
**Endolite India Ltd.**