Echelon^{VAC}

User Guide

EVAC22L1S—EVAC30R8S
EVAC22L1SD—EVAC30R8SD

EN	User Guide	2
SL	Uporabniški priročnik	9
BG	Наръчник на потребителя	16
HR	Korisnički priručnik	23
SK	Používateľská príručka	30
HU	Felhasználói útmutató	37
EL	Οδηγός χρήσης	44
LV	Lietotāja rokasgrāmata	51
LT	Naudotojo vadovas	58
ET	Kasutaja käsiraamat	65

- Contents2
- 1 Description and Intended Purpose3
- 2 Safety Information3
- 3 Construction4
- 4 Maintenance6
- 5 Limitations on Use6
- 6 Technical Data7

1 Description and Intended Purpose

Please keep these instructions.

The term device is used throughout this document to refer to EchelonVAC .

Application

This device is to be used exclusively as part of a lower limb prosthesis.

Intended for single user.

The device provides limited self-alignment of the prosthesis on varied terrain and following footwear changes. In addition to the visco-elastic, self aligning hydraulic ankle it generates a vacuum in the range 12-17 in Hg. It is intended to improve on postural sway and symmetry while easing abnormal pressures at the socket interface. A moderate energy return foot with resilient multi-axial ankle movement. Independent heel and toe spring provides some axial deflection and the split toe provides good ground compliance.

2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Any changes in the performance or function of the limb e.g. restricted movement, non-smooth motion or unusual noises should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

The device is not suitable for extreme sports, running or cycle racing, ice and snow sports, extreme slopes and steps. Any such activities undertaken are done so completely at the users' own risk. Recreational cycling is acceptable.

Assembly, maintenance and repair of the device must only be carried out by a suitably qualified clinician.

The user should contact their clinician if their condition changes.

Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.

To minimise the risk of slipping and tripping, appropriate footwear that fits securely onto the footshell must be used at all times.

After continuous use the ankle casing may become hot to the touch.

Avoid exposure to extreme heat and/ or cold.

The user must not adjust or tamper with the setup of the device.


Be aware of finger trap hazard at all times.

3 Construction


Principal Parts:

- Hydraulic Body Assembly including pyramid (Aluminium/St. Stl./Titanium)
- Carrier Assembly (Aluminium/St. Stl.)
- Heel & Toe Springs (e-Carbon)
- Spring Attachment Screws (Titanium/St.Stl.)
- Glide Sock (UHM PE)
- Foot Shell (PU)
- Vacuum Parts (PU, Nylon, Aluminium)



Vacuum System Parts

- 1a Auto-Expulsion Valve
- 1b Threaded Housing
- 2 Check Valve
- 3 In-line Filter
- 4 Vacuum Tubing

Vacuum System

Vacuum system guidance

	<i>Symptoms:</i>	<i>Cause/Remedy</i>
1.	Unable to generate a vacuum	Vacuum tube(s) disconnected Inspect and reconnect as necessary Limited ankle movement due to footwear a full range of movement is required
2.	Unable to maintain a vacuum	Vacuum tube(s) disconnected Inspect and reconnect as necessary Leakage at socket valve/barbs Consult your practitioner Check the integrity of the vacuum seal at the socket/residuum interface

4 Maintenance

Maintenance must be carried out by competent personnel.

We recommend a regular visual check at least annually. Signs of wear that may affect function should be reported to their service provider (e.g. significant wear or excessive discoloration from long term exposure to UV).

Any changes in performance of this device must be reported to the practitioner.

Changes in performance may include:

- Increase in ankle stiffness
- Reduced ankle support (free movement)
- Any unusual noise
- Lack of vacuum

Note: Check that the suspension sleeve maintains a good seal with the skin proximally and socket distally.

The practitioner must also be informed of:

- Any changes in body weight and/or activity level.
- Discoloration of the residual limb.
- Significant wear or excessive discoloration from long term exposure to UV.

Cleaning

Use a damp cloth and mild soap to clean outside surfaces. Do not use aggressive cleansers.

5 Limitations on Use

Intended life:

A local risk assessment should be carried out based upon activity and usage.

Lifting loads:

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

Environment:

Avoid exposing the device to corrosive elements such as water, acids and other liquids. Also avoid abrasive environments such as those containing sand for example, as these may promote premature wear.

Exclusively for use between -15°C and 50°C (5°F to 122°F).


Suitable for outdoor use

6 Technical Data

Operating and Storage Temperature Range:	-15 °C to 50 °C (5 °F to 122 °F)
Component Weight [Size 26N]:	930 g (2 lb 1 oz)
Recommended Activity Level:	2, 3, 4
Maximum User Weight:	125 kg (275 lb)
Proximal Alignment Attachment:	Male Pyramid (Blatchford)
Range of Hydraulic Ankle Motion: (excludes additional range of motion provided by heel and toe springs)	6 degrees plantar flexion to 3 degrees dorsiflexion
Build Height:	[Sizes 22-24] 120 mm
[See diagram below]	[Sizes 25-26] 125 mm
	[Sizes 27-30] 130 mm
Heel Height:	10 mm
Maximum Vacuum:	17in Hg

Fitting length

Size	A
22-24	120 mm
25-26	125 mm
27-30	130 mm

Size	B
22-26	65 mm
27-28	70 mm
29-30	75 mm

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: www.blatchford.co.uk

Medical Device

Single Patient – multiple use

Warranty

This device is warranted for 36 months - foot shell 12 months - glide sock 3 months.

This warranty does not apply to:

Consumable parts including the vacuum tubing, in-line filter and valves unless a failure has occurred due to a defect in materials or workmanship.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

See Blatchford website for the current full warranty statement.

Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

Environmental Aspects

Where possible the components should be recycled in accordance with local waste handling regulations.

Trademark Acknowledgements

Echelon and Blatchford are registered trademarks of Blatchford Products Limited.

Manufacturer's Registered Address

Blatchford Products Limited, Lister Road, Basingstoke RG22 4AH.

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© Blatchford Products Limited 2021. All rights reserved.
Endolite India Ltd.
Blatchford Europe GmbH