AqualimbTF

Instructions for Use

AQTF22L-AQTF27R

EN	Instructions for Use	2
DA	Brugsanvisning	15
NO	Bruksanvisning	28
FI	Käyttöohjeet	41
SV	Bruksanvisning	54
RU	Инструкция протезиста	67
ZH	使用说明	80
AR	إرشادات الاستخدام	93

Contents	2
1 Description and Intended Purpose	3
2 Safety Information	4
3 Construction	5
4 Function.....	6
5 Maintenance	7
6 Limitations on Use	8
7 Assembly Instructions	9
8 Fitting Instructions.....	10
9 Bench Alignment.....	11
10 Dynamic Alignment	12
11 Technical Data	12
12 Ordering Information	13

1 Description and Intended Purpose

These Instructions for Use are intended for use by the practitioner and user unless otherwise stated.

The term *device* is used throughout these instructions for use to refer to AqualimbTF.

Please read and ensure you understand all instructions for use, in particular all safety information and maintenance instructions.

Application

This device is to be used exclusively as part of a lower limb prosthesis.

Intended for a single user.

This device is a prosthetic knee designed for above-knee users.

This device is water-resistant and is intended for use **only** as a shower limb.

Features

- Geometrically stable 4-bar knee
- Hand operated knee lock (HOKL)
- Corrosion resistance

Activity Level

This device is for Activity Level 1 and 2 users. Activity Level 3 and 4 users may only use this device if they limit their activities to Activity Level 2.

Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low-level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

Contraindications

- This device should not be used for any activity above Activity Level 2.
- Users with insufficient voluntary control should only use the device locked.
- Not for use by transtibial users.

Clinical Benefits

- Allows the user to wear their prosthesis whilst showering and in wet environments.

2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Any changes in performance (for example 'play', instability, lock not engaging), or changes in swing resistance should immediately be reported to your service provider.

Take all necessary safety precautions and due care on wet and slippery surfaces.

The user should be advised to contact their practitioner if their condition changes.

The user must not adjust or tamper with the setup of the device.

Avoid exposure to extreme heat and/or cold.

The device is not suitable for extreme sports, running or cycle racing, ice and snow sports, extreme slopes and steps. Any such activities undertaken are done so completely at the users' own risk.

Be aware of finger trap hazard at all times

Assembly, maintenance and repair of the device must only be carried out by a suitably qualified practitioner.

The freedom of swing of the knee will vary slightly between when the knee is dry and when it is wet and /or contaminated with soap for example.


The device is designed for prolonged submersion and suitable for immersion in fresh water only. Ensure any use of the device in water complies with the conditions given in Section 6 *Limitations on Use*.

3 Construction

Principal Parts

- Knee Mechanism Acetal homopolymer, stainless steel
- Top Housing Glass fiber reinforced nylon
- Foot Keel/Shin Glass fiber reinforced nylon/PU
- Attachment Parts Stainless steel/ aluminum alloy

Component Identification

4 Function

The knee has a polycentric design comprising pivots connecting a chassis to a main body via three link arms. A hand-operated spring-loaded pin operates within a latched slot in the body. The pin may be located to engage into the chassis and lock the knee in full extension. A housing is bolted to the distal aspect of the knee into which is clamped a foot/shin component. The shin may be shortened to suit the user. The sole of the foot features an anti-slip tread pattern.

Knee Lock Operation


If lock mode is required ensure that the lock actuator is fully engaged before walking or showering.

The lock is hand operated. To lock the knee release the actuator from the down position, fully extend the knee and push the actuator up to the top position and ensure the actuator is fully engaged as shown.

Knee Unlocked

Actuator in down position


Knee shown in free mode.

The lock will not engage at full knee extension.

Knee Locked

Actuator fully engaged in top position

Knee shown with lock in place.

5 Maintenance

Visually check the device regularly.

Report any changes in performance of this device to the practitioner/service provider e.g. unusual noises, increased stiffness or restricted/excessive rotation, significant wear or excessive discoloration from long term exposure to UV.

Inform the practitioner/service provider of any changes in body weight and/or activity level.

Cleaning

Use a damp cloth and mild soap to clean outside surfaces. **Do not** use aggressive cleansers, wipe dry.

The remaining instructions in this section are for practitioner use only.

This maintenance must be carried out only by competent personnel (practitioner or suitable trained technician).

The following routine maintenance is to be carried out at least annually:

- Check the security of the proximal attachment nut and distal alignment bolt.
- Check for satisfactory function of knee lock and freedom of swing.
- Visually check for wear. Some surface damage may occur after a period of use; this does not affect the function or strength of the device.

The user should be advised:

- Any changes in performance of this device must be reported to the practitioner.

Changes in performance may include:

- 'Play' in moving parts
- Lock not engaging
- Instability

Check condition of knee bumper and replace if necessary.

Ensure the user has read and understood all safety and user-level maintenance information.

Advise the user that a regular visual check of the device is recommended and signs of wear that may affect function (e.g. significant wear or excessive discoloration from long term exposure to UV) and/or any changes in performance should be reported to their service provider.

Advise the user to inform the practitioner/service provider of any changes in body weight and/or activity level.

6 Limitations on Use

Intended Life

Service life of the product is covered by the warranty period — a local risk assessment should be carried out based upon activity and usage.

Lifting Loads

User weight and activity is governed by the stated limits. Load carrying by the user should be based on a local risk assessment.

Environment

This device is waterproof up to a maximum depth of 1 meter.

To prevent damage or wear, thoroughly rinse this device with fresh water after use in the following environments:

- Abrasive environments, such as those that contain sand or grit;
- Salt or chlorinated water.

Finish this device to prevent water ingress (see *Assembly Instructions*). If water enters it, invert it to allow the water to drain, then leave it to dry.

Exclusively for use between -15 °C and 50 °C (5 °F and 122 °F).

Suitable for submersion

7 Assembly Instructions

The instructions in this section are for practitioner use only.

-  Use appropriate health and safety equipment at all times including extraction facilities.
-  Be aware of finger trap hazard at all times.

1

Manufacture the socket using the appropriate Blatchford single bolt former. Cut a 27 mm diameter hole in distal end of socket. Place socket plate inside socket. Assemble alignment parts supplied. Minimum socket thickness 6 mm.

2

Apply Loctite 243 and torque the proximal alignment bolt to 35 Nm. Use a peg spanner (940060) to prevent the socket plate from rotating.

8
35 Nm

3

Apply Loctite 243 to bolt threads. Assemble the distal top housing using the nut and washer as shown and torque tighten to 35 Nm after alignment.


4

It may be easier to fit the shin before assembling the top housing. Ensure full engagement of shin, see diagram Section *Bench Alignment*.

8 Fitting Instructions


The instructions in this section are for practitioner use only.

1


Fit shin to housing ensuring full engagement. Measure length to determine amount to be cut from shin.

2


Determine final cut length required.

3


Remove shin from housing, cut to length and reassemble then check measurement.

4

Tighten clamp screws evenly, applying a torque of 5 Nm at least twice to each screw. Example order of tightening shown to evenly spread load.


5

Check length and alignment. Finalise the screw settings during dynamic alignment.

9 Bench Alignment

The instructions in this section are for practitioner use only.

10 Dynamic Alignment


Set up as shown in *Bench Alignment* if the user requires adjustment length, shift, rotation and angular adjustments are available by loosening, adjusting and re-tightening of distal and proximal alignment fixings.

Ensure after any adjustments are made, that the load line is not posterior to the top pivot pin.

11 Technical Data

Operating and Storage Temperature Range:	-15 °C to 50 °C 5 °F to 122 °F
Total Component Weight (<i>size 24</i>):	2 kg (4.4lb)
Activity Level:	1-2
Maximum User Weight:	100 kg
Proximal Alignment Attachment:	Blatchford single bolt alignment
Build Height: (See diagram below)	400–585 mm

Fitting Length

12 Ordering Information

Order Example

AQ	TF	25	L
	Transfemoral	Size	Side (L - left/R - right)

Available from size 22 to size 27:

AQTF22L–AQTF27R

e.g. AQTF25L, AQTF26R

Replacement Parts

Part	Part No.
Aqualimb TF Proximal Alignment Kit	169205
Knee Bumper	622310
Screw Proximal M10 x 55 ST STL	910420
Peg Spanner	940060

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: www.blatchford.co.uk

Medical Device

Single Patient – multiple use

Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a practitioner.

Warranty

This device is warranted for 24 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

See the Blatchford website for the current full warranty statement.

Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

Environmental Aspects

Where possible the components should be recycled in accordance with local waste handling regulations.

Retaining the Packaging Label

The practitioner is advised to keep the packaging label as a record of the device supplied.

Trademark Acknowledgements

Blatchford is a registered trademark of Blatchford Products Limited.

Manufacturer's Registered Address

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Blatchford Europe GmbH
Endolite India Ltd.