AqualimbTT

Instructions for Use

629705L-629716L

EN	Instructions for Use	2
DA	Brugsanvisning	14
NO	Bruksanvisning	26
FI	Käyttöohjeet	38
SV	Bruksanvisning	50
RU	Инструкция протезиста	62
ZH	使用说明	74
AR	إرشادات الاستخدام	86

Contents

Contents	2
1 Description and Intended Purpose	3
2 Safety Information	4
3 Construction	5
4 Function.....	6
5 Maintenance	6
6 Limitations on Use	7
7 Bench Alignment.....	8
7.1 Static Alignment	8
7.2 Dynamic Alignment	8
8 Assembly Instructions	9
9 Technical Data	11
10 Ordering Information	12

1 Description and Intended Purpose

These Instructions for Use are intended for use by the practitioner and user unless otherwise stated.

The term *device* is used throughout these instructions for use to refer to AqualimbTT.

Please read and ensure you understand all instructions for use, in particular all safety information and maintenance instructions.

Application

This device is to be used exclusively as part of a lower limb prosthesis.

Intended for a single user.

This device is a shin-foot unit that is for use as a component of a transtibial prosthesis.

This device is water-resistant and is intended for use only as a shower limb.

Activity Level

This device is for Activity Level 1 and 2 users. Activity Level 3 and 4 users may only use this device if they limit their activities to Activity Level 2.

Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low-level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

Contraindications

This device should not be used for any activity above Activity Level 2.

The device is not to be used with footwear.

Users with poor balance (unless they remain seated while showering).

Clinical Benefits

Allows the user to wear their prosthesis whilst showering and in wet environments.

2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Any changes in performance should immediately be reported to your service provider.

Take all necessary safety precautions and due care on wet and slippery surfaces.

The user should be advised to contact their practitioner if their condition changes.

The user must not adjust or tamper with the setup of the device.

Avoid exposure to extreme heat and/or cold.

The device is not suitable for extreme sports, running or cycle racing, ice and snow sports, extreme slopes and steps. Any such activities undertaken are done so completely at the users' own risk.

Be aware of finger trap hazard at all times

Assembly, maintenance and repair of the device must only be carried out by a suitably qualified practitioner.

This device is buoyant. If the user intends to submerge this device, remove the cosmesis with a suitable cutting tool (see *Assembly Instructions*).


The device is designed for prolonged submersion and suitable for immersion in fresh water only. Ensure any use of the device in water complies with the conditions given in *Limitations on Use*.

3 Construction

Principal Parts

- Cosmesis PU foam
- Foot Shell PU foam
- Shin Glass-fiber-reinforced nylon
- Alignment Kit Stainless steel/aluminum/titanium/
glass-fiber-reinforced nylon
- Chassis Glass-fiber-reinforced nylon
- M5x50 Socket Cap Screw Stainless steel
- Clamp Plates Stainless steel

Component Identification

4 Function

This device is a shin-foot unit that comprises a foot shell with an integrated shin. The shin can be shortened to suit the user's height. The sole of the foot features an anti-slip tread for grip on wet surfaces.

5 Maintenance

Visually check the device regularly.

Report any changes in performance of this device to the practitioner/service provider e.g. unusual noises, increased stiffness or restricted/excessive rotation, significant wear or excessive discoloration from long term exposure to UV.

Inform the practitioner/service provider of any changes in body weight and/or activity level.

Cleaning

Use a damp cloth and mild soap to clean outside surfaces. **Do not** use aggressive cleansers, wipe dry.

The remaining instructions in this section are for practitioner use only.

The following routine maintenance is to be carried out at least annually:

- Make sure that all screws are tight. If they are loose, tighten them to the correct torque setting;
- Make sure that the anti-slip tread functions as intended, and replace if necessary.

Advise the user to report the following changes to the practitioner:

- Changes in body weight;
- Changes in the performance of the device, for example, unusual noises or play.

Ensure the user has read and understood all safety and user-level maintenance information.

Advise the user that a regular visual check of the device is recommended and signs of wear that may affect function (e.g. significant wear or excessive discoloration from long term exposure to UV) and/or any changes in performance should be reported to their service provider.

Advise the user to inform the practitioner/service provider of any changes in body weight and/or activity level.

6 Limitations on Use

Intended Life

A local risk assessment should be carried out based upon activity and usage.

Lifting Loads

User weight and activity is governed by the stated limits. Load carrying by the user should be based on a local risk assessment.

Environment

This device is waterproof up to a maximum depth of 1 meter.

To prevent damage or wear, thoroughly rinse this device with fresh water after use in the following environments:

- Abrasive environments, such as those that contain sand or grit;
- Salt or chlorinated water.

Finish this device to prevent water ingress (see *Assembly Instructions*). If water enters it, invert it to allow the water to drain, then leave it to dry.

Exclusively for use between -15 °C and 50 °C (5 °F and 122 °F).

Suitable for submersion

This device is buoyant. If the user intends to submerge this device, remove the cosmetics with a suitable cutting tool (see *Assembly Instructions*).

Take care on wet and slippery surfaces.

7 Bench Alignment

The instructions in this section are for practitioner use only.

7.1 Static Alignment

Sagittal Plane

* With the flexion accommodated, the load line should pass through the target shown in the image, ± 10 mm A-P. The target intersects the point where the anterior face of the shin meets the keel.

Make sure that when the user is standing, the heel and toe are evenly loaded and the foot is in full contact with the floor.


7.2 Dynamic Alignment

Coronal Plane

Minimize M-L thrust by adjusting the relative positions of the socket and the foot .

Sagittal Plane

Adjust the relative position of the socket and the foot to ensure that the user transitions smoothly from heel strike to toe-off.







8 Assembly Instructions

The instructions in this section are for practitioner use only.

 Use appropriate health and safety equipment at all times including extraction facilities.

 Be aware of finger trap hazard at all times.

8 Assembly Instructions (cont.)

The instructions in this section are for practitioner use only.


7

Go to *Bench Alignment* section, then:
continue to step 8 or 9.

8




If the user intends to submerge this device, remove the cosmesis with a suitable cutting tool.

Note: This action is irreversible.

9


If the user does not intend to submerge this device, finish it to prevent water ingress.

9 Technical Data

Operating and Storage Temperature Range:	-15 °C to 50 °C (5 °F to 122 °F)
Component Weight (size 26):	1.39 kg (3 lb)
Activity Level:	1-2
Maximum User Weight:	100 kg (220 lb)
Build Height:	See diagram below

Fitting Length

10 Ordering Information

Order Example

Size (cm)	Part Number — Left	Part Number — Right
22	629705L	629706L
23	629707L	629708L
24	629709L	629710L
25	629711L	629712L
26	629713L	629714L
27	629715L	629716L

Replacement Parts

Part	Part Number
Alignment Kit	169106
M5x50 Socket Cap Screw	910415

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: www.blatchford.co.uk

Medical Device

Single Patient – multiple use

Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a practitioner.

Warranty

This device is warranted for 24 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

See the Blatchford website for the current full warranty statement.

Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

Environmental Aspects

Where possible the components should be recycled in accordance with local waste handling regulations.

Retaining the Packaging Label

The practitioner is advised to keep the packaging label as a record of the device supplied.

Trademark Acknowledgements

Blatchford is a registered trademark of Blatchford Products Limited.

Manufacturer's Registered Address

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Blatchford Europe GmbH
Endolite India Ltd.