Dynamic Pylon Top Housing

Instructions for Use

539601
539602

| | | |
|----|----------------------------|----|
| EN | Instructions for Use | 2 |
| FR | Instructions d'utilisation | 11 |
| DE | Gebrauchsanweisung | 20 |
| IT | Istruzioni per l'uso | 29 |
| ES | Instrucciones de uso | 38 |
| NL | Gebruiksaanwijzing | 47 |
| PL | Instrukcje użytkowania | 56 |
| PT | Instruções de utilização | 65 |
| CS | Návod k použití | 74 |

Contents

| | |
|--|---|
| Contents | 2 |
| 1 Description and Intended Purpose | 3 |
| 2 Safety Information | 4 |
| 3 Construction | 5 |
| 4 Function..... | 6 |
| 5 Maintenance | 6 |
| 6 Limitations on Use | 7 |
| 7 Fitting Advice..... | 8 |
| 8 Technical Data | 9 |
| 9 Ordering Information | 9 |

1 Description and Intended Purpose

These Instructions for Use are intended for the practitioner and user unless otherwise stated.

The term *device* is used throughout to refer to the Dynamic Pylon Top Housing.

Please read and ensure you understand these instructions, in particular all safety information and maintenance instructions.

Application

This device is to be used exclusively for modular lower limb prostheses, intended for a single user.

The device is designed to connect to other Blatchford-approved prosthetic technology (Elite Blade or Javelin), providing secure connection with sliding and rotational adjustment.

Features

- Lightweight, strong, aluminum and titanium construction

Activity Level

This device is suitable for Activity Levels 1-4 (weight limits apply, see *Technical Data*). Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances and any such decision should be made with sound and thorough justification.

Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low-level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

Contraindications

There are no known contraindications if used in accordance with these instructions.

Clinical Benefit

- Allows a connection from a blade to other Blatchford-approved prosthetic technology.

2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Any changes in the performance or function of the limb e.g. restricted movement, non-smooth motion, excessive play or unusual noises should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

539601 (male) is designed for prolonged submersion and suitable for immersion in fresh water only. Ensure any use of the device in water complies with the conditions given in *Limitations on Use*.

539602 (female) can be used as shower prosthesis but is not designed for prolonged submersion in water. After contact with water wipe dry immediately. If this device comes into contact with salt or chlorinated water, it should be rinsed with fresh water and dried. Ensure any use of the device in water complies with the conditions given in *Limitations on Use*.

The user must not adjust or tamper with the setup of the device.

The user is advised to contact their practitioner if their condition changes.

To reduce the risk of injury due to failure or loosening of the bolt connections, ensure the bolt threads are cleaned thoroughly before each installation.

Always use Loctite and apply the specified torque value to the bolts. Never use an alternative bolt.


Assembly, maintenance and repair of the device must only be carried out by a suitably qualified practitioner.

3 Construction

Principal Parts

- Female Pyramid Titanium
- Pyramid Set Screws Plated steel (x4)
- Collar/Clamp Ring Aluminum alloy
- Clamp Screw and Inserts Plated steel
- Male Pyramid Stainless steel
- Dome Aluminum
- Pylon Top Housing Assembly Aluminum
- Fixing Bolts Titanium (x2)

Component Identification

4 Function

This device provides secure connection between a male or female pyramid and a transtibial blade prosthesis. When the fixing bolts are tightened securely, using Loctite 243, this device offers a secure way, with sliding, rotating and angular movement, to connect and optimize alignment of the two parts of the limb.

5 Maintenance

Visually check the device regularly.

Report any changes in performance of this device to the practitioner/service provider e.g. unusual noises or significant wear. There should be no movement at this interface.

Inform the practitioner/service provider of any changes in body weight and/or activity level.

Cleaning

Use a damp cloth and mild soap to clean outside surfaces. DO NOT use aggressive cleansers.

ALWAYS dry the device thoroughly after shower use to prevent water collection in the “T” slot of the Pylon Top Housing Assembly.

The remaining instructions in this section are for practitioner use only.

This maintenance must be carried out only by competent personnel (practitioner or suitable trained technician).

The following routine maintenance is to be carried out at least annually:

- Check the fixing bolts for tightness. If loose, remove and clean, then apply Loctite 243 and tighten to the correct torque setting, see *Construction* section.
- Check for defects that could affect proper function.
- Check for corrosion.

Ensure the user has read and understood all safety and user-level maintenance information.

Advise the user that a regular visual check of the device is recommended and signs of wear that may affect function should be reported to their service provider (e.g. significant wear or corrosion).

Advise the user to inform the practitioner/service provider of any changes in body weight and/or activity level.

If this device is used for extreme activity, the maintenance level and interval should be reviewed and if required advice and technical support sought to plan a new maintenance schedule dependent upon the frequency and nature of the activity. This should be determined by a local risk assessment carried out by a suitably qualified individual.

6 Limitations on Use

Intended Life

A local risk assessment should be carried out based upon activity and usage.

Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

Environment

539601 (Male)

This device is waterproof to a maximum depth of 1 meter. Thoroughly rinse with fresh water after use in abrasive environments such as those that may contain sand or grit, for example, to prevent wear or damage to moving parts. Thoroughly rinse with fresh water after use in salt or chlorinated water.

ALWAYS dry the device thoroughly after shower use to prevent water collection in the "T" slot of the Pylon Top Housing Assembly.

Exclusively for use between -15 °C and 50 °C (5 °F to 122 °F).

Suitable for submersion

539602 (Female)

This device can be used as a shower prosthesis but is not designed for prolonged submersion in water. If this device comes into contact with salt or chlorinated water, it should be rinsed with fresh water and dried. Light surface corrosion will not affect the function or security of this device. However, if heavy corrosion is evident, stop using this device and contact your practitioner.

ALWAYS dry the device thoroughly after shower use to prevent water collection in the "T" slot of the Pylon Top Housing Assembly.

Exclusively for use between -15 °C and 50 °C (5 °F to 122 °F).

Suitable for shower use

7 Fitting Advice


The instructions in this section are for practitioner use only.

A fitting jig (941255) is available to optimize length and alignment. Alignment of the Prosthesis can then be finalized for comfort, function and performance in conjunction with the required Male/Female pyramid adapter. Position as close as possible to the build line.

Shift Adapter

Refer to appropriate Dynamic Pylon fitting instructions for cutting/drilling information.

Insert relevant shift adapter into 'T' slot, if required. Attach housing, ensuring pylon does not protrude above the top surface of the adapter.

Fixing Bolts

Two fixing bolts are supplied with the kit for fixing to mating component.


Apply Loctite 243 to all screws threads and tighten to the correct torque setting.

| Symptom | Solution |
|---|---|
| A recurring noise occurs between the device interfaces. | Clean the device interfaces. Apply Loctite and tighten the bolts to the correct torque setting. |
| The adapter moves out of position. | User must not use the device until adjusted, repaired or replaced. |

8 Technical Data

| | 539601 | 539602 |
|--|---|---|
| Operating and Storage Temperature Range: | -15°C to 50°C (5°F to 122°F) | -15°C to 50°C (5°F to 122°F) |
| Component Weight: | 165 g (5.8 oz) | 240 g (8.5 oz) |
| Activity Level: | 1–4 | 1–4 |
| Maximum User Weight: | Levels 1–3 :166 kg (366 lb) Level 4: 130 kg (287 lb) | Levels 1–3 :166 kg (366 lb) Level 4: 130 kg (287 lb) |
| Range of Adjustment: | ±7° tilt from vertical ±15 mm shift A-P | ±7° tilt from vertical ±15 mm shift A-P |
| Proximal Attachment: | Male pyramid (Blatchford) | Female pyramid (Blatchford) |
| Build Height: | 7 mm | 28 mm |

Fitting Length

Storage and Handling

When storing for prolonged periods, ensure the product is free from moisture and stored at room temperature.

9 Ordering Information

| Device | Part Number |
|--|-------------|
| Dynamic Pylon Top Housing (Male) Round Profile | 539601 |
| Dynamic Pylon Top Housing (Female) Round Profile | 539602 |

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: www.blatchford.co.uk

Medical Device

Single Patient – multiple use

Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

Warranty

This device is warranted for 24 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

See the Blatchford website for the current full warranty statement.

Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

Environmental Aspects

This product is made from a recyclable material. Where possible, the components should be recycled in accordance with local waste handling regulations.

Retaining the Packaging Label

The practitioner is advised to keep the packaging label as a record of the device supplied.

Trademark Acknowledgements

Blatchford is a registered trademark of Blatchford Products Limited.

Manufacturer's Registered Address

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Blatchford Europe GmbH
Endolite India Ltd.