Proximal Pyramid Adapter Kit

Instructions for Use

239017, 239089

EN	Instructions for Use	2
SL	Navodila za uporabo	15
BG	Инструкции за употреба	28
HR	Upute za uporabu	41
SK	Návod na použitie	54
HU	Használati útmutató	67
EL	Οδηγίες χρήσης	80
LV	Lietošanas pamācība	93
LT	Naudojimo instrukcija	106
ET	Kasutusjuhend	119

Contents

Contents	2
1 Description and Intended Purpose	3
2 Safety Information	4
3 Construction	5
4 Function.....	6
5 Maintenance	6
6 Limitations on Use	7
7 Bench Alignment.....	7
7.1 Static Alignment.....	7
7.2 Dynamic Alignment	7
8 Fitting Advice.....	8
8.1 Dynamic Pylon Products	8
8.2 Mercury Knee	9
8.3 ESK+ 4 Bolt Chassis.....	10
8.4 Linx, Orion3	10
8.5 KX06V2/ KX07	11
9 Technical Data	12
10 Ordering Information	12

1 Description and Intended Purpose

These Instructions for Use are intended for use by the practitioner and user unless otherwise stated.

The term *device* is used throughout these instructions for use to refer to the Proximal Pyramid Adapter Kit.

Please read and ensure you understand all instructions for use, in particular all safety information and maintenance instructions.

Application

The device is designed for use exclusively as part of a lower limb prosthesis.

The device is approved for use with Blatchford modular components.

Intended for single user only.

Activity Level

This device is suitable for activity levels 1 to 4; weight limits apply, see *Technical Data*, Section 8. Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances and any such decision should be made with sound and thorough justification.

Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited ambulator.

Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

Activity Level 3

Has the ability or potential for ambulation with variable cadence.

Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

Contraindications

- Only for use with Blatchford devices.
- 239089 - Not for use if rotation required.
- Use only with the dome in the correct orientation.

Clinical Benefits

- **Allows a connection between the prosthetic knee or pylon to other prosthetic componentry, which can be customised to optimise alignment.**

2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Any changes in the performance or function of the device e.g. unusual noises should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

Avoid exposure to extreme heat and/or cold.

To reduce the risk of injury due to failure or loosening of the screw connections ensure the bolt threads are cleaned thoroughly before each installation.

The device is designed for prolonged submersion and suitable for immersion in fresh water only. Ensure any use of the device in water complies with the conditions given in Section 6 *Limitations on Use*.

Be aware of finger trap hazard at all times.

Assembly, maintenance and repair of the device must only be carried out by a suitably qualified clinician that has attended an approved training course.

The user must not adjust or tamper with the setup of the device.

The user should be advised to contact their practitioner if their condition changes.

Always apply the specified torque value to bolt. Never replace the bolt with an alternative bolt. Always use the bolt specified.


Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.

3 Construction

Principal Parts

- Pyramid
 - Dome
 - O-ring
- Titanium
Aluminum
Nitrile Rubber

239017

239089

4 Function

Allows attachment between a Blatchford pyramid system and Blatchford components with “T” slot.

Allows A-P slide and rotation when used in conjunction with a Blatchford “T” slot.

The non-rotating version of the device allows for slide function only.

5 Maintenance

Visually check the device regularly.

Report any changes in performance of this device to the practitioner/service provider.

Changes in performance may include:

- Any unusual noise
- Any play or loss alignment

Inform the practitioner/service provider of any changes in body weight and/or activity level.

Cleaning

Use a damp cloth and mild soap to clean outside surfaces, do not use aggressive cleansers.

The remaining instructions in this section are for practitioner use only.

Maintenance must be carried out only by competent personnel (practitioner or suitable trained technician).

Carry out the following routine maintenance at least annually:

- Ensure all screws are secure. If not remove and clean screws, reapply Loctite and tighten all screws to the correct torque settings.
- Check for visual defects that may affect proper function.

Ensure the user has read and understood all safety and user-level maintenance information.

Advise the user that a regular visual check of the device is recommended and signs of wear, deterioration or corrosion that may affect function should be reported to their service provider.

Advise the user to inform the practitioner/service provider of any changes in body weight and/or activity level.

6 Limitations on Use

Intended Life

A local risk assessment should be carried out based upon activity and usage.

Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

Environment

This device is waterproof to a maximum depth of 1 meter.

Thoroughly rinse with fresh water after use in abrasive environments such as those that may contain sand or grit, for example, to prevent wear or damage to moving parts.

Thoroughly rinse with fresh water after use in salt or chlorinated water.

Suitable for submersion

Exclusively for use between -15 °C and 50 °C
(5 °F to 122 °F).

7 Bench Alignment

The instructions in this section are for practitioner use only.

7.1 Static Alignment

Adjust the alignment (A-P slide and tilt (and rotation 239017)) in combination with interfacing adapters to properly accommodate flexion, adduction or abduction. Ensure that the wearer can stand in a stable position with the heel and toe are evenly loaded.

7.2 Dynamic Alignment

Coronal Plane

Ensure that M-L thrust is minimal by adjusting relative positions of socket and foot.

Sagittal Plane

Adjust relative positions of Prosthetic components to achieve a smooth, stable transition from heel strike to toe-off.

Note... Ensure that the wearer can stand in a stable position with the heel and toe are evenly loaded.

8 Fitting Advice

The instructions in this section are for practitioner use only.

Avoid extremes of alignment.

To ensure the correct range of A-P slide ensure the dome is fitted the correct way around or orientated correctly. Ensure the anti-rotation feature of 239089 is assembled correctly.

To avoid loss of alignment in use, ensure interfacing pyramid alignment screws are torqued and secured correctly.

Position as close to the build line as possible, if possible avoid extremes of alignment.

Key to following Diagrams


- * = Anti-rotation feature posterior
- ** = Anti-rotation feature anterior

8.1 Dynamic Pylon Products


239017

239089



938306




8.2 Mercury Knee

239017


938251

239089



8.3 ESK+ 4 Bolt Chassis

239017


938205

239089

8.4 Linx, Orion3


239017

Linx

938353


239089

Orion3



938374




8.5 KX06V2/ KX07

239017

938409


239089

9 Technical Data

Operating and Storage Temperature Range:	-15°C to 50°C (5°F to 122°F)
Component Weight:	50 g
Activity Level:	1–4
Maximum User Weight:	166 kg (366 lb)
Build Height:	7 mm

Fitting Length

10 Ordering Information

Part Numbers

239017 Male Pyramid Adapter–Rotating

239089 Male Pyramid Adapter–Non-rotating

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: www.blatchford.co.uk.

Medical Device

Single Patient – multiple use

Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

Warranty

This device is warranted for 24 months.

Foot shell is warranted for 12 months.

Glide sock is warranted for 3 months.

This warranty does not apply to consumable parts including the vacuum tubing, in-line filter and valves unless a failure has occurred due to a defect in materials or workmanship.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

See the Blatchford website for the current full warranty statement.

Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

Environmental Aspects

This product is fabricated from a polyurethane that cannot be easily recycled: please dispose of it responsibly as general waste, as per local waste handling regulations.

Where possible the components should be recycled in accordance with local waste handling regulations.

Retaining the Packaging Label

You are advised to keep the packaging label as a record of the device supplied.

Trademark Acknowledgements

Avalon and Blatchford are registered trademarks of Blatchford Products Limited.

Manufacturer's Registered Address

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Blatchford Europe GmbH
Endolite India Ltd.