# TT Pro

## Instructions for Use

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EN	Instructions for Use	2
DA	Brugsanvisning	16
NO	Bruksanvisning	30
FI	Käyttöohjeet	44
SV	Bruksanvisning	58
RU	Инструкция протезиста	72
ZH	使用说明	86
AR	إرشادات الاستخدام	100

# Contents

Contents .....	2
1 Description and Intended Purpose .....	3
2 Safety Information .....	4
3 Construction .....	5
4 Function.....	6
5 Maintenance .....	6
6 Limitations on Use .....	6
7 Bench Alignment.....	7
8 Fitting Advice.....	8
9 Dis-assembling TT Pro Pyramid.....	8
10 Assembling TT Pro Pyramid.....	9
11 Dis-assembling TT Pro Clamp .....	10
12 Assembling TT Pro Clamp .....	11
13 Technical Data .....	12
14 Ordering Information .....	14

---

# 1 Description and Intended Purpose

These Instructions for Use are intended for use by the practitioner and user unless otherwise stated.

The term device is used throughout these instructions for use to refer to TT Pro.

Please read and ensure you understand all instructions for use, in particular all safety information and maintenance instructions.

## Application

This device is for use as a component of a lower limb prosthesis.

Intended for single user.

The device is a shock and torsion device that balances shock absorption and energy management. Reduced stress at the interface between the body and the prosthesis provides greater comfort for users.

## Activity Level

This device is recommended for users that have the potential to achieve Activity Level 3 or 4. Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances and any such decision should be made with sound and thorough justification.

Note... The maximum user weight for Level 4 users is 100 kg (220 lb).

### Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

### Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low-level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

### Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

### Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

## Contraindications

This device may not be suitable for Activity Level 1 individuals or for competitive sports events, as these types of users will be better served by a specially designed prosthesis optimized for their needs.

Ensure that the user has understood all the instructions for use, drawing particular attention to the maintenance and safety information.

## Clinical Benefits

- Reduced loading rate on the residual limb, compared to rigid pylons
- Reduced back pain when twisting, compared to rigid pylons
- Increased walking speed, compared to rigid pylons


# Torsion Rod and Spring Selection

## Torsion Rod

1 = 1 dot

2 = 2 dots

3 = 3 dots

## Spring

P = Purple

W = White

B = Black

Weight kg	Impact Level		
	Low	Moderate	High
44–52	1P	2P	3W
53–59	2P	2W	3W
60–68	2P	3W	3B
69–77	2W	3W	3B
78–88	2W	3B	3B
89–100	3B	3B	3B
101–116	3B	3B	—
117–125	3B	3B	—

## 2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Any changes in the performance or function of the limb e.g. restricted movement, non-smooth motion or unusual noises should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

The device is not suitable for extreme sports, running or cycle racing, ice and snow sports, extreme slopes and steps. Any such activities undertaken are done so completely at the users' own risk. Recreational cycling is acceptable.

Assembly, maintenance and repair of the device must only be carried out by a suitably qualified clinician.

The user should contact their clinician if their condition changes.

Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.

To minimize the risk of slipping and tripping, appropriate footwear that fits securely onto the footshell must be used at all times.


# 3 Construction

## Principal Parts


- Top Housing/ Tube Clamp Assembly (aluminum, carbon steel, stainless steel, acetal homopolymer, PTFE, PU)
- Torsion Rod (acetal homopolymer)

- Shaft/Locking Ring Assembly (stainless steel/aluminum alloy)
- Spring (carbon steel)
- Shin Assembly (stainless steel/aluminum)

### TT Pro Pyramid

### TT Pro Clamp

---

## 4 Function

The device comprises an energy efficient coil spring and a polymer torsion rod. Together, these components reduce both shear force at the socket interface during heel strike and torsional loads during turning. Energy that is stored in the axial spring during stance phase is returned at toe-off to give a smooth energy efficient gait. The bearing surfaces are coated with titanium nitride, which provides a hard-wearing low friction surface for moving parts.

---

## 5 Maintenance

The user should be advised to report the following to the practitioner:

- Changes in body weight or activity level
- Changes in the performance of this device, for example:
  - Unusual noises or play.
  - Reduced energy return.

If this device is used for extreme activity, the maintenance level and interval should be reviewed and if required advice and technical support sought to plan a new maintenance schedule dependent upon the frequency and nature of the activity. This should be determined by a local risk assessment carried out by a suitably qualified individual.

### Cleaning

Use a damp cloth and mild soap to clean outside surfaces. **Do not** use aggressive cleansers.

*The remaining instructions in this section are for practitioner use only.*

Only competent personnel should carry out maintenance.

Carry out the following maintenance annually:

- Disassemble, clean, and lubricate parts as shown in dis-assembly and assembly instruction Sections 9-12.
- Check for rotational play
- Replace the torsion rod if it's worn

---

## 6 Limitations on Use

### Intended Life

A local risk assessment should be carried out based upon activity and usage.

### Lifting Loads

User weight and activity is governed by the stated limits. Load carrying by the user should be based on a local risk assessment.

### Environment

Avoid exposing this device to corrosive elements such as water, acids and other liquids. Avoid abrasive environments such as those containing sand, as these may promote premature wear.


Exclusively for use between -15 °C and 50 °C (5°F and 122 °F).

Suitable for outdoor use

## 7 Bench Alignment

*The instructions in this section are for practitioner use only.*

### Additional Components

Follow the alignment instructions that are supplied with any additional components.

The build line should pass through the center of the device.

\* Allow for the user's footwear.

## 8 Fitting Advice

*The instructions in this section are for practitioner use only.*

The device uses an independent axial spring and Torsion rod to control compressive and torsional resilience. Typical axial movement should be between 5-8 mm for normal walking.

When used in combination with other devices that have an axial compression that give an overall shortening effect on the limb, the overall deflection may be excessive. In these cases, a stiffer axial spring may be appropriate. If the axial spring is changed, function should be checked to confirm that torsional compliance has been maintained.

Range of rotational movement can be affected by alignment. Excessive torsional movement coupled with instability in late stance is indicative of insufficient M-L stability and M-L positioning of the foot in relation to the weight line. If instability is lateral, reposition the foot medially, conversely if instability is medially, reposition the foot laterally. These adjustments should be checked before considering a stiffer torsion rod.

## 9 Dis-assembling TT Pro Pyramid




*The instructions in this section are for practitioner use only.*

**Be aware of finger trap hazard at all times.**






**Use appropriate health and safety equipment at all times including extraction facilities.**





# 10 Assembling TT Pro Pyramid

The instructions in this section are for practitioner use only.







# 11 Dis-assembling TT Pro Clamp

*The instructions in this section are for practitioner use only.*

# 12 Assembling TT Pro Clamp

The instructions in this section are for practitioner use only.

---



## 13 Technical Data

Operating and Storage Temperature Range:	-15 °C to 50 °C (5 °F to 122 °F)
Component Weight:	Pyramid— 523 g (18 oz) Clamp— 440 g (15 oz)
Activity Level:	1–4
Maximum User Weight:	Activity Level 1–3— 125 kg (275 lb) Activity Level 4— 100 kg (220 lb)
Range of Adjustment:	Proximal Attachment: $\pm 7^\circ$ Angular Distal Attachment: 360 ° Rotation
Build Height:	TT Pro Pyramid— 135–226 mm TT Pro Clamp— 109 mm
Proximal Attachment:	Female Pyramid (Blatchford)
Distal Attachment:	TT Pro Pyramid— Blatchford $\varnothing 30$ mm Tube Clamp TT Pro Clamp— Blatchford $\varnothing 30$ mm Tube


### Fitting Length: TT Pro Pyramid

Note... Fit the device with the screws facing anteriorly.

### Fitting Length: TT Pro Clamp

Note... Fit the device with the screws facing anteriorly.

# 14 Ordering Information

## Order Example


### TT PRO Pyramid

TTPRO-PYR3B

### TT PRO Clamp

TTPRO-30-3B

Part	Part Number
Torsion Rod	406007 0.5Nm (1 dot)
	406008 0.6Nm (2 dots)
	406009 0.7Nm (3 dots)
Spring	406014 Purple (P)
	406015 White (W)
	406016 Black (B)

Part	Part Number
M8x16 Flat Point Socket Set Screw	910474
M8x10 Flat Point Socket Set Screw	910476
M3x10 Socket Cap Screw	910012
M4x4 Cup Point Socket Set Screw	910013
M6x20 Socket Cap Screw	910048
M6 Lock Washer	916150
Spring Cap Washer	406013
Thrust Washer	405031
Progressive Stop Washer	406004
30 mm Shin Tube Cutting Jig	941237
TT Locking Ring Spanner	940070

## Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for damage caused by the component combinations that were not authorized by the manufacturer.

## CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: [www.blatchford.co.uk](http://www.blatchford.co.uk)

Medical Device

Single Patient – multiple use

## Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

## Warranty

This device is warranted for 24 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

See the Blatchford website for the current full warranty statement.

## Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

## Environmental Aspects

This product is constructed from recyclable materials. Where possible the components should be recycled in accordance with local waste handling regulations.

## Retaining the Packaging Label

The practitioner is advised to keep the packaging label as a record of the device supplied.

## Manufacturer's Registered Address

Blatchford Products Limited, Lister Road, Basingstoke RG22 4AH, UK.

## Trademark Acknowledgements

TT Pro and Blatchford are registered trademarks of Blatchford Products Limited.

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**Blatchford Europe GmbH**
**Endolite India Ltd.**