Mercury

Instructions for Use

MH30 / MH30 SWING (*as shown*)

MHPYR / MHPYR SWING

MHDPYRS, MHDPYRM, MHD30S, MHD30M

239452

| | | |
|----|-----------------------|----|
| EN | Instructions for Use | 2 |
| DA | Brugsanvisning | 12 |
| NO | Bruksanvisning | 22 |
| FI | Käyttöohjeet | 32 |
| SV | Bruksanvisning | 42 |
| RU | Инструкция протезиста | 52 |
| ZH | 使用说明 | 62 |
| AR | إرشادات الاستخدام | 72 |

Contents

- Contents2
- 1 Description and Intended Purpose3
- 2 Safety Information5
- 3 Construction6
- 4 Function.....7
- 5 Maintenance7
- 6 Limitations on Use8
- 7 Alignment and Fitting Advice.....8
- 8 Technical Data9
- 9 Ordering Information 10

1 Description and Intended Purpose

Application

The term *device* is used throughout these instructions for use to refer to Mercury knee products. These instructions are for the practitioner.

This device is a knee unit that provides a uniaxial frame for Swing and Stance position activated support for the prosthesis for walking at various cadence on varied terrain (a swing only version is also available). It is intended to improve safety and reduce compensation.

Ensure that the user understands all instructions for use, drawing particular attention to the maintenance and safety information.

Features

| | Swing & Stance | Swing Only |
|--|----------------|------------|
| Tough carbon fibre frame | ✓ | ✓ |
| 125° knee flexion | ✓ | ✓ |
| Long life needle roller bearings | ✓ | ✓ |
| Reinforced hydraulic unit | ✓ | ✓ |
| Durable urethane kneeler pad | ✓ | ✓ |
| Variable cadence | ✓ | ✓ |
| Compact design | ✓ | ✓ |
| Adjustable stance yield resistance | ✓ | |
| Swing only mode (suitable for cycling) | ✓ | ✓ |
| Cylinder flexion lock | ✓ | |

Clinical Benefits

Swing and Stance (S&S)

- Flexion / Extension adjustability allows walking at variable speeds.
- Knee resistance yield allows supported stair descent.
- Knee resistance yield allows supported sitting.
- Reduced shock at heel strike compared to a locking knee.
- Flexion lock allows extended standing support.
- Selectable swing only mode to suit some sports activities'.

Swing Only

- Flexion / Extension adjustability allows walking at variable speeds.
- Maximum level of voluntary control for competent users.

Activity Level

This device is for Activity Level 3 and 4 users.

Contraindications

This device might not be suitable for Activity Level 1–2 users or for use in competitive sports events. These types of users might be better served by a specially designed prosthesis that is optimized for their needs.

This device is not suitable for users with:

1. Residual muscular weakness, contractures that cannot be correctly accommodated, or proprioceptive dysfunction including poor balance.
2. Contralateral joint instabilities or pathology.
3. Complicated conditions involving multiple disabilities.

This device is for a single user.

Users must be given gait training before using this device.

2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Be aware of finger trap hazard at all times.

Any changes in the performance or function of the limb e.g. instability, double-action, restricted movement, non-smooth motion or unusual noises should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

Any excessive changes in heel height after finalization of alignment may adversely affect limb function.

After continuous use the cylinder can become hot.

Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.

Care should be taken when carrying heavy loads as this may adversely affect the stability of the device.

The user should be advised to contact their practitioner if their condition changes.

Avoid exposure to extreme heat and/or cold.


Assembly, maintenance and repair of the device must only be carried out by a suitably qualified clinician.

3 Construction

Principal Parts

- Carrier Assembly - CF. Composite & St/Steel
- Chassis Assembly - Alu Alloy & St/Steel
- Hydraulic Control Cylinder - Various inc. Alu Alloy, St/Steel & Hydraulic Oil
- Distal Dome - Alu Alloy
- Distal Pyramid - Titanium with St/Steel bolt
- Distal Tube Clamp - Alu Alloy
- Knee Pad - Thermoplastic Elastomer

For assembly of Cylinder to Carrier Assembly 239452 & 239037/38/40/41 please see document 938218.

4 Function

The knee is a uniaxial design comprising a knee pivot connecting a chassis to a carrier. A hydraulic knee control sits within this frame assembly, fitted between the chassis posterior pivot and the carrier distal pivot.

The S&S device has an adjustable yielding hydraulic support activated from knee extension prior to heel strike and released by a hyper-extension moment at any time, usually occurring after mid-stance.

The swing only device has no stance control mechanism and requires a combination of geometric stability and voluntary control during stance phase.

Both devices provide an adjustable hydraulic swing phase control to accommodate variable cadence.

5 Maintenance

An annual visual inspection is recommended. Check for visual defects that may affect proper function.

Maintenance must be carried out by competent personnel.

The user should be advised:

- Any changes in performance of this device must be reported to the practitioner.

Changes in performance may include:

- Increase or decrease in knee stiffness.
- Instability.
- Reduced knee support (free movement).
- Any unusual noise.*

* Due to the nature of the hydraulics there may be some slight air noise from the cylinder during the first few steps. This is not detrimental to the function of the unit and should dissipate quickly. If symptoms persist, please consult your practitioner.

Storage and Handling

When storing for prolonged periods place the device vertically with trunnion uppermost. Alternatively and in addition, air management may be assisted by cycling the knee control several times with the stance resistance switched off.

Use product packaging supplied.

Cleaning

Use a damp cloth and mild soap to clean the outside surfaces.

DO NOT use aggressive cleansers.

6 Limitations on Use

Intended Life

A local risk assessment should be carried out based upon activity and usage.

Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

Environment

Avoid exposing the device to corrosive elements such as water, acids and other liquids. Also avoid abrasive environments such as those containing sand for example as these may promote premature wear.

Suitable for outdoor use

Exclusively for use between -10 °C and 50 °C (14 °F to 122 °F).


It is recommended that only Blatchford products be used in conjunction with the device.

7 Alignment and Fitting Advice

For alignment and fitting advice please see document 938216.


8 Technical Data

| | |
|--|--|
| Operating and Storage Temperature Range: | -10 °C to 50 °C (14 °F to 122 °F) |
| Component Weight: | 1.23Kg (2lb 11oz) |
| Activity Level: | 3, 4 |
| Maximum User Weight: | 150Kg (330lb) - Activity L3 145Kg (319lb) - Activity L4 |
| Maximum Flexion Angle | 125 ° |
| Proximal Alignment attachment: | 4-bolt interface T-Slot Interface |
| Distal Alignment attachment: | Blatchford Male Pyramid OR 30mm Shin Tube Clamp |
| Build Height | See Diagram |

For applications above 125kg (275lb) the appropriate shin kit MUST be used.

*See Section 9

9 Ordering Information

| Description | Part Number |
|--|-------------|
| <i>Device</i> | |
| Mercury Hi-Activity Knee PYR | MHPYR |
| Mercury Hi-Activity Knee PYR Small K/Ball dis-con | MHDPYRS |
| Mercury Hi-Activity Knee PYR Medium K/Ball dis-con | MHDPYRM |
| Mercury Hi-Activity 30mm tube clamp | MH30 |
| Mercury Hi-Activity 30mm tube vclamp Small K/Ball dis-con | MHD30S |
| Mercury Hi-Activity 30mm tube clamp Medium K/Ball dis-con | MHD30M |
| Mercury Hi-Activity Knee Carrier Assy less cyl/distal attachment | 239452 |
| <i>Accessories *</i> | |
| High strength shin kit double pyramid | 239029 |
| High strength shin kit single pyramid | 239030 |

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: www.blatchford.co.uk

Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

Medical Device

Single Patient – multiple use

Warranty

This device is warranted for 36 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

See the Blatchford website for the current full warranty statement.

Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

Environmental Aspects

This product contains hydraulic oil, mixed metals and plastics. It should be recycled where possible in accordance with local waste recycling regulations.

Retaining the Packaging Label

You are advised to keep the packaging label as a record of the device supplied.

Trademark Acknowledgements

Blatchford is a registered trademark of Blatchford Products Limited.

Manufacturer's Registered Address

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Blatchford Europe GmbH
Endolite India Ltd.