# Navigator

## Instructions for Use

NAV22LWW-NAV30RBB  
NAV22LWWD-NAV30RBBD

|    |                            |     |
|----|----------------------------|-----|
| EN | Instructions for Use       | 2   |
| FR | Instructions d'utilisation | 16  |
| DE | Gebrauchsanweisung         | 30  |
| IT | Istruzioni per l'uso       | 44  |
| ES | Instrucciones de uso       | 58  |
| NL | Gebruiksaanwijzing         | 72  |
| PL | Instrukcje użytkowania     | 86  |
| PT | Instruções de utilização   | 100 |
| CS | Návod k použití            | 114 |

# Contents

|  |    |
|--|----|
| Contents .....                                   | 2  |
| 1 Description and Intended Purpose .....         | 3  |
| 1.1 Ankle Ball and Snubber Selection .....       | 4  |
| 2 Safety Information .....                       | 5  |
| 3 Construction .....                             | 6  |
| 4 Function.....                                  | 7  |
| 5 Maintenance.....                               | 7  |
| 6 Limitations on Use .....                       | 8  |
| 7 Alignment .....                                | 9  |
| 7.1 Additional Components .....                  | 9  |
| 7.2 Static Alignment.....                        | 9  |
| 7.3 Dynamic Alignment .....                      | 9  |
| 8 Fitting Advice.....                            | 10 |
| 9 Assembly Instructions .....                    | 10 |
| 9.1 Changing the Ankle Ball .....                | 10 |
| 9.2 Changing the Snubber.....                    | 11 |
| 9.3 Assembling and Finishing the Foot Shell..... | 12 |
| 10 Technical Data .....                          | 13 |
| 11 Ordering Information .....                    | 14 |

---

# 1 Description and Intended Purpose

These instructions are for the practitioner and user unless otherwise stated.

The term *device* is used throughout these instructions for use to refer to Navigator.

Please read and ensure you understand all instructions for use, in particular all safety information and maintenance instructions.

## Application

This device is an ankle-foot unit that is for use as a component of a lower limb prosthesis.

This device comprises an integrated multi-axial ankle and a resilient keel which offer the user enhanced ground compliance.

This device is for a single user.

## Features

- Multi-axial ankle with resilient keel
- Short, energy-efficient, rollover shape keel optimized for K2 gait
- Integrated ankle joint
- Choice of ankle stiffness
- Sandal toe foot shell

## Activity Level

This device is recommended for Activity Level 1–3 users (weight restrictions apply, see *Technical Data* section). As an exception to our recommendation, Activity Level 4\* users may only use this device based upon the outcome of a prosthetic assessment.

### Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

### Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low level environmental barriers such as, curbs, stairs or uneven surfaces. Typical of the limited community ambulator.

### Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

### Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

\*Maximum user weight: 100 kg (220 lb).

## Contraindications

This device is not suitable for use in competitive sports events.

## Clinical Benefits

- Optimized keel allows for more consistent rollover radius of curvature.
- Rollover shape falls within the most energy efficient range for most people.
- Greater stability for low mobility users when moving through stance phase compared to single axis feet and solid-ankle-cushioned-heel feet.
- Improved inter-limb symmetry of stance phase timing, compared to rigid ankles.
- Reduced sound limb loading, compared to rigid ankles.

For bilateral users, the clinical benefits of this device compared to rigid ankle devices include:

- No loss of stability during standing.
- Easier to walk on uneven ground, on slopes and to walk longer distances.

## 1.1 Ankle Ball and Snubber Selection


The ankle ball and the snubber are color-coded to indicate their hardness.

Use the table to choose an ankle ball and snubber combination that suits the weight and Activity Level of the user. The letter on the left refers to the ankle ball. For example, WO refers to a white ankle ball and an orange snubber.

- White (W) = Soft
- Orange (O) = Medium
- Blue (B) = Hard

Notes:

- Bilateral users might need stiffer snubbers than unilateral users.
- Transfemoral users might need softer ankle balls than transtibial users.

| Weight  |         | Activity Level |    |    |
|---------|---------|----------------|----|----|
| kg      | lb      | 1              | 2  | 3  |
| 44–52   | 100–115 | WW             | WW | WW |
| 53–59   | 116–130 | WW             | WW | WO |
| 60–68   | 131–150 | WW             | WO | OO |
| 69–77   | 151–170 | WO             | OO | OB |
| 78–88   | 171–195 | OO             | OB | BB |
| 89–100  | 196–220 | OB             | BB | BB |
| 101–125 | 221–275 | BB             | BB | BB |

## 2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Any changes in the performance or function of the limb e.g. restricted movement, non-smooth motion or unusual noises should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

The device is not suitable for extreme sports, running or cycle racing, ice and snow sports, extreme slopes and steps. Any such activities undertaken are done so completely at the users' own risk. Recreational cycling is acceptable.

Assembly, maintenance and repair of the device must only be carried out by a suitably qualified clinician.

The device is not intended for use when immersed in water or as a shower prosthesis. If the limb comes into contact with water wipe dry immediately. Ensure any use of the device complies with the conditions given in *Limitations on Use*.

Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.

To minimize the risk of slipping and tripping, appropriate footwear that fits securely onto the foot shell must be used at all times.

Avoid exposure to extreme heat and/or cold.

The user must not adjust or tamper with the setup of the device.

The user should contact their clinician if their condition changes.


Be aware of finger trap hazard at all times.

# 3 Construction

## Principal Parts

- Dome (aluminum alloy)
- Snubber (rubber)
- Top Housing (aluminum alloy)
- Ankle Ball and Pyramid (aluminum alloy/rubber/stainless steel/titanium)
- Keel (fiber-reinforced polymer)
- M6x27 Socket Cap Screw (plated carbon steel)
- Glide Sock (UHM PE)
- Foot Shell (polyurethane)

## Component Identification

---

## 4 Function

This device is an ankle-foot unit that comprises both an integrated ankle and a resilient keel inside a sandal-toe foot shell. The ankle center is positioned anatomically to promote a natural gait from heel strike to toe-off. The plantar flexion and dorsiflexion resistances can be adjusted by changing the ankle ball and the snubber to suit the weight and the Activity Level of the user.

---

## 5 Maintenance

Visually check the device regularly.

Report any changes in performance of this device to the practitioner/service provider e.g. unusual noises, increased stiffness or restricted/excessive rotation, significant wear or excessive discoloration from long term exposure to UV.

Inform the practitioner/service provider of any changes in body weight and/or activity level.

### **Cleaning**

Use a damp cloth and mild soap to clean outside surfaces, DO NOT use aggressive cleansers.

***The remaining instructions in this section are for practitioner use only.***

This maintenance must be carried out only by competent personnel (practitioner or suitable trained technician).

The following routine maintenance is to be carried out at least annually:

- Check the foot shell and glide sock for damage or wear, and replace if necessary;
- Check the ankle ball and the snubber for damage or wear, and replace if necessary;
- Make sure the socket cap screws are tight. If they are loose, tighten them to the correct torque setting (see *Construction* section);
- Re-lubricate the snubber.

Ensure the user has read and understood all safety and user-level maintenance information.

Advise the user that a regular visual check of the device is recommended and any changes in performance and signs of wear that may affect function should be reported to their service provider (e.g. significant wear or excessive discoloration from long term exposure to UV).

Advise the user to report to the practitioner any changes in body weight or activity level.

---

## 6 Limitations on Use

### Intended Life

A local risk assessment should be carried out based upon activity and usage.

### Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

### Environment

Avoid exposing the device to corrosive elements such as water, acids and other liquids. Avoid abrasive environments such as those containing sand, as these may promote premature wear.


Suitable for outdoor use

Exclusively for use between -15 °C and 50 °C  
(5 °F and 122 °F).

## 7 Alignment

### 7.1 Additional Components

Follow the alignment instructions that are supplied with any additional components.

### 7.2 Static Alignment

#### Sagittal Plane

With the flexion accommodated, the load line should pass through the center of the socket and through the center of the pyramid.

Make sure that, when the user is standing, the heel and the toe are evenly loaded and the foot is in full contact with the floor.

### 7.3 Dynamic Alignment

#### Coronal Plane

Minimize M-L thrust by adjusting the relative positions of the socket and the foot.

#### Sagittal Plane

Make sure that the user transitions smoothly from heel strike to toe-off.

\* Accommodate the heel height of the user's footwear

## 8 Fitting Advice

The main function of the ankle ball is to provide plantar flexion resistance. Harder ankle balls provide more resistance and softer ankle balls provide less resistance.

The main function of the snubber is to provide dorsiflexion resistance. Harder snubbers provide more resistance and softer snubbers provide less resistance.

| Symptoms  | Remedy  |
|---|---|
| <ul style="list-style-type: none"><li>• Sinking at heel strike</li><li>• Foot slap</li></ul>  | <ol style="list-style-type: none"><li>a) Confirm A–P alignment and heel height accommodation.</li><li>b) Fit a harder ankle ball.</li></ol> |
| <ul style="list-style-type: none"><li>• Rapid transition from heel strike through stance phase</li><li>• Difficulty controlling heel action (e.g. the foot jars into mid-stance)</li><li>• Foot feels too rigid</li></ul> | <ol style="list-style-type: none"><li>a) Confirm A–P alignment and heel height accommodation.</li><li>b) Fit a softer ankle ball.</li></ol> |
| <ul style="list-style-type: none"><li>• Drop-off</li><li>• Feels like walking downhill</li></ul>  | <ol style="list-style-type: none"><li>a) Confirm A–P alignment and heel height accommodation.</li><li>b) Fit a harder snubber.</li></ol>    |
| <ul style="list-style-type: none"><li>• Feels like walking uphill</li></ul>   | <ol style="list-style-type: none"><li>a) Confirm A–P alignment and heel height accommodation.</li><li>b) Fit a softer snubber.</li></ol>    |


Please contact your supplier if it is not possible to achieve a smooth gait after following the advice above.

## 9 Assembly Instructions


 **Be aware of finger trap hazard at all times.**

 **Use appropriate health and safety equipment at all times including extraction facilities.**

### 9.1 Changing the Ankle Ball

1. To release ankle ball, remove keel and ankle assembly from foot shell and remove four M6x27 socket cap screws (see Fig. 3).
2. Clean the replacement ankle ball (Fig. 1) and the keel pocket with acetone or a similar solvent.
3. Align the flats of the male pyramid along the axis of the keel and rest the ankle ball in the keel pocket.
4. Insert the alignment screw into the ankle ball through the hole in the bottom of the keel.


5. Clean inside of the top housing with acetone or a similar solvent.

6. Place the top housing over the ankle ball.

Note... Make sure the opening of the top housing faces posteriorly.

7. Insert four socket cap screws fully into the four holes in the bottom of the keel, then apply Loctite 243 to the screw threads where they expose above the keel.

8. Apply pressure to the top housing while maintaining pyramid alignment. Then, engage the socket cap screws into the housing.


9. Turn each screw, in sequence, ½ to 1 turn at a time until the top housing comes into contact with keel.

Note... Make sure that the alignment screw stays in place while you tighten the socket cap screws.

10. Torque tighten the socket cap screws to 10 Nm and remove the alignment screw.

## 9.2 Changing the Snubber


1. Lubricate the snubber with a suitable rubber lubricant.

2. Press the cone tool onto the housing.

3. Place the dome over the ankle stem.

Note... Make sure the 'A' on top of the dome faces anteriorly.

## 9.3 Assembling and Finishing the Foot Shell

1. Place the keel inside the glide sock and slide them both into the foot shell.

Note... Make sure that the heel of the keel is fully engaged in the location feature in the heel of the foot shell.

Note... After assembling to the limb, make sure that the glide sock is not trapped under the dome.

2. Pull the glide sock under slight tension.
3. Use a cable tie to secure the glide sock to the pylon just above the adapter.

4. Remove the end of the cable tie, then fold the glide sock over.


5. Abrade the top surface of the foot shell.
6. Bond a 12 mm pelite interface to the foot shell and shape as appropriate. Use Evostik 528 (926213) or an equivalent.
7. Bond a foam cosmesis to the top surface of the pelite interface and shape as appropriate. Use Thixofix (926204), Evostik 528 (926213) or an equivalent.

# 10 Technical Data

|  |                                     |
|--|-------------------------------------|
| Operating and Storage Temperature Range: | -15 °C to 50 °C<br>(5 °F to 122 °F) |
| Component Weight (size 26):              | 565 g (20 oz)                       |
| Activity Level:                          | 1-3                                 |
| Maximum User Weight:                     | 125 kg (275 lb)                     |
| Proximal Alignment Attachment:           | Male pyramid (Blatchford)           |
| Range of Adjustment                      | ±7° angular                         |
| Build Height:                            | 95 mm                               |
| Heel Height:                             | 10 mm                               |

## Fitting Length

# 11 Ordering Information

## Order Example

|            |           |               |                          |                    |               |
|------------|-----------|---------------|--------------------------|--------------------|---------------|
| <b>NAV</b> | <b>25</b> | <b>L</b>      | <b>W</b>                 | <b>W</b>           | <b>D</b>      |
|            | Size      | Side<br>(L/R) | Ankle<br>Ball<br>(W/O/B) | Snubber<br>(W/O/B) | *Dark<br>Tone |

Available from size 22 to size 30:

NAV22LWW to NAVR30RBB

NAV22LWWD to NAV30RBBD

\*(Add 'D' for a dark tone foot shell)

e.g. NAV25LOB, NAV29RWW, NAV26RBBD

### Ankle and Keel Assemblies (USA only)

| Size        | Type       | Part Number |
|-------------|------------|-------------|
| Small       | White 60°  | 409156      |
| Medium      | White 60°  | 409157      |
| Large       | White 60°  | 409158      |
| Extra Large | White 60°  | 409159      |
| Small       | Orange 70° | 409256      |
| Medium      | Orange 70° | 409257      |
| Large       | Orange 70° | 409258      |
| Extra Large | Orange 70° | 409259      |
| Small       | Blue 80°   | 409356      |
| Medium      | Blue 80°   | 409357      |
| Large       | Blue 80°   | 409358      |
| Extra Large | Blue 80°   | 409359      |

### Snubbers

| Type       | Part Number |
|------------|-------------|
| White 60°  | 405605      |
| Orange 70° | 405606      |
| Blue 80°   | 405607      |

Glide Sock 405815

### Foot Shells (Add a D for dark tone)

| Size (cm) | Left Part Number | Right Part Number |
|-----------|------------------|-------------------|
| 22        | 539020           | 539021            |
| 23        | 539022           | 539023            |
| 24        | 539024           | 539025            |
| 25        | 539026           | 539027            |
| 26        | 539028           | 539029            |
| 27        | 539030           | 539031            |
| 28        | 539032           | 539033            |
| 29        | 539034           | 539035            |
| 30        | 539036           | 539037            |

### Tools for Assembly

| Part            | Part Number |
|-----------------|-------------|
| Alignment Screw | 910414      |
| Cone Tool       | 940083      |
| Cable Tie       | 929029      |

### Ankle Ball and Pyramids

| Type       | Part Number |
|------------|-------------|
| White 60°  | 409456      |
| Orange 70° | 409556      |
| Blue 80°   | 409656      |

## Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

## CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: [www.blatchford.co.uk](http://www.blatchford.co.uk)

Medical Device

Single Patient – multiple use

## Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

## Warranty

This device is warranted for 24 months, the foot shell is warranted for 12 months, and the glide sock is warranted for 3 months. The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions. See Blatchford website for the current full warranty statement.

## Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

## Environmental Aspects

Where possible the components should be recycled in accordance with local waste handling regulations.

## Retaining the Packaging Label

You are advised to keep the packaging label as a record of the device supplied.

## Trademark Acknowledgements

Navigator and Blatchford are registered trademarks of Blatchford Products Limited.

## Manufacturer's Registered Address

Blatchford Products Limited, Lister Road, Basingstoke RG22 4AH, UK

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**Blatchford Europe GmbH**
**Endolite India Ltd.**