SmartIP

Instructions for Use

SMARTIP4B30
SMARTIP4BPYR

SMARTIPSF30
SMARTIPSFPYR

EN	Instructions for Use	2
SL	Navodila za uporabo	24
BG	Инструкции за употреба	46
HR	Upute za uporabu	68
SK	Návod na použitie	90
HU	Használati útmutató	112
EL	Οδηγίες χρήσης	134
LV	Lietošanas pamācība	156
LT	Naudojimo instrukcija	178
ET	Kasutusjuhend	200

Contents	2
1 Description and Intended Purpose	3
2 Safety Information	4
3 Construction	5
4 Function.....	6
5 Maintenance	6
6 Limitations on Use	7
7 Static Alignment	8
8 Limb Setup	9
8.1 SmartIP Operation	9
8.2 Fine Tuning.....	9
8.3 Cushion Valve Adjustment.....	10
8.4 Using the Reset Switch.....	11
9 Alignment Advice and Stability Adjustment.....	12
10 Fitting Advice.....	13
10.1 Changing Stanceflex Rubbers	13
10.2 Front Stop Adjustment - Adjustable Trunnion	14
10.3 Further Adjustments.....	15
11 Assembly Instructions	16
11.1 Changing Front Stop Buffer.....	16
11.2 Changing Shin Stops.....	16
11.3 Removing ESK+ Mechanism	17
11.4 Replacement of ESK+ Mechanism	18
12 Technical Data	19
13 Ordering Information	20

1 Description and Intended Purpose

These Instructions for Use are intended for use by the practitioner and user unless otherwise stated.

The term *device* is used throughout these instructions for use to refer to SmartIP.

Please read and ensure you understand all instructions for use, in particular all safety information and maintenance instructions.

Application

This device is to be used exclusively as part of a lower limb prosthesis.

Intended for a single user.

The device is a programmable microprocessor controlled knee with adjustable terminal impact, extension damping and weight activated stance control.

Features

- Smart Programming for flexion control
- Manually adjustable terminal impact and extension control
- Comfortable gait from very slow to fast walking speeds
- Smooth speed progression because of enhanced program
- Weight activated stance control
- Stanceflex option improves comfort
- Battery life expectancy 9-12 months under normal usage
- 140° knee flexion
- 4-Bolt option reduces build height

Activity Level

This device is recommended for users that have the potential to achieve Activity Level 2 or 3. Some Activity Level 4 users may also benefit from using the device (weight limits apply to all activity levels, see *Technical Data*).

Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances and any such decision should be made with sound and thorough justification.

Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low-level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

Activity Level 3

Has the ability or potential for ambulation with variable cadence.

Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

Contraindications

This device may not be suitable for Activity Level 1.

Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

Clinical Benefits

Microprocessor controlled swing device (IP) proven to:

- Increase stability and mobility on slopes and stairs, compared to non microprocessor-controlled knees
- Reduce energy expenditure and cognitive demand, compared to non microprocessor-controlled knees
- Make walking easier at different speeds, compared to non microprocessor-controlled knees

2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Beware of finger trap hazard all times.

Any changes in the performance or function of the limb e.g. instability, double-action, restricted movement, non-smooth motion or unusual noises should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

Any excessive changes in the heel height after finalization of alignment may adversely affect limb function.

Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.

Care should be taken when carrying heavy loads as this may adversely affect the stability of the device.

The user should be advised to contact their practitioner if their condition changes.

Avoid exposure to extreme heat and/or cold.

Assembly, maintenance and repair of the device must only be carried out by a suitably qualified clinician.

The user must not adjust or tamper with the setup of the device.

The device is not intended for use when immersed in water or as a shower prosthesis. If the limb comes into contact with water wipe dry immediately.

Avoid strong magnetic fields, sources of electrical interference, atmospheres containing liquids and/or powders.


Contact your service provider for advice if the device is to be used with a pacemaker or any other electronic medical device.

3 Construction

Principal Parts

- Carrier Assembly Carbon fibre composite & stainless steel
- Chassis Assembly Aluminium alloy, stainless steel, engineering thermoplastic, fibre reinforced (asbestos free) friction material.
- Distal Dome Aluminium alloy
- Distal Pyramid Titanium with stainless steel bolt
- Distal Tube Clamp Aluminium alloy
- Patella Engineering thermoplastic
- Stanceflex Rubber Natural rubber
- Stanceflex strap Kevlar
- SmartIP Unit Aluminium alloy, stainless steel, Titanium, engineering thermoplastic Electronic components, batteries

Component Identification

4 Function

The device is a uniaxial design comprising a knee pivot connecting a chassis to a carrier. A weight activated stabilizing knee mechanism is housed within the chassis. The stabilizing mechanism has a resilient action which provides limited stance flexion under load to absorb shock loads at heel strike.

The Stanceflex knee chassis has a pivotally connected top plate which is resiliently supported by a rubber block, this provides an increased stance flexion action under load to further absorb shock loads at heel strike.

A knee control cylinder sits within this frame assembly, fitted between the chassis posterior pivot and the carrier distal pivot.

The SmartIP control cylinder uses a sensors to monitor walking speed and automatically adjusts flexion resistance accordingly. Extension control is via fixed adjustment of a terminal impact damper.

5 Maintenance

Visually check the device regularly.

Report any changes in performance of this device to the practitioner/service provider e.g. unusual noises, increased stiffness, restricted/excessive rotation or significant wear.

Inform the practitioner/service provider of any changes in body weight and/or activity level.

Cleaning

Use a damp cloth and mild soap to clean outside surfaces, do not use aggressive cleansers.

The remaining instructions in this section are for practitioner use only.

This maintenance must be carried out only by competent personnel (practitioner or suitable trained technician).

To maintain optimum performance, an annual inspection for unacceptable noise, play and stability is recommended. More frequent inspections may be required depending on factors such as the activity of the user, frequency of usage and the stability setting of the knee.

Ensure the user has read and understood all safety and user-level maintenance information.

Advise the user that a regular visual check of the device is recommended and any changes in performance of this device must be reported to the practitioner. Changes in performance may include: increase or decrease in knee stiffness, instability, reduced knee support (free movement), any unusual noise.

If this device is used for extreme activity, the maintenance level and interval should be reviewed and if required advice and technical support sought to plan a new maintenance schedule dependent upon the frequency and nature of the activity. This should be determined by a local risk assessment carried out by a suitably qualified individual.

6 Limitations on Use

Intended Life

A local risk assessment should be carried out based upon activity and usage.

Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

Environment

Avoid exposing the device to corrosive elements such as water, acids and other liquids. Also avoid abrasive environments such as those containing sand for example as these may promote premature wear.

Exclusively for use between -10 °C to 50 °C
(14 °F to 122 °F)


Suitable for outdoor use

7 Static Alignment


The instructions in this section are for practitioner use only.

1. Complete static alignment as per diagram below.

2. Walking Trials.

⚠ ESX+ function must be confirmed by wearer before walking on the limb. (See Alignment Advice and Stability Adjustment)

Fully extend limb to access the adjuster.

Rotate the SPRING STACK ADJUSTER to achieve the required stability for the individual amputee.

Note: Audible clicks should be heard upon adjustment, and the adjuster is self locking.

Tighten the SPRING STACK to reduce stability.

Release the SPRING STACK to increase stability - See above.


If required, the adjuster may be reset to the factory setting. This can be achieved by reducing stability by six audible clicks from the maximum stability setting.

8 Limb Setup

Adjust for static alignment as shown in the ESK+ fitting instructions supplied with the module.
Cushion Valve adjustment - If after programming for the first time, the user experiences too much or too little damping at the end of the swing, the cushion valve may be adjusted. (See *Cushion Valve Adjustment*) Fine Tuning MUST always follow any cushion valve adjustment. (See *Fine Tuning*)


8.1 SmartIP Operation

The instructions in this section are for practitioner use only.

8.2 Fine Tuning

The instructions in this section are for both practitioner and user.

8.3 Cushion Valve Adjustment

The instructions in this section are for practitioner use only.

Adjust the Cushion Valve prior to commencement of programming.


For increased damping of terminal impact, rotate the valve clockwise (2.5 A/F Wrench Key).

Adjust the cushion resistance whilst the amputee is walking slowly.

At this speed the damping level should be set to the highest resistance that still permits the knee to extend fully and remain extended during the entire stance phase. This should result in the elimination of terminal impact at higher speeds.

Fine tune settings as the trial progresses, always confirming that the limb fully extends during slow walking.

Note... Typical increments are 1/4 turn until valve has effect thereafter smaller increments should be used to optimize terminal impact.

8.4 Using the Reset Switch

The instructions in this section are for practitioner use only.

Initial Set Up, Calibration and Programming

- Selects factory default program settings
 - Enters calibration mode followed by programming mode.
1. Press reset button, confirmed by 1 beep.
 2. Within 2 seconds press reset button again, confirmed by 2 beeps.
 3. Within 2 seconds press reset button again, confirmed by 3 beeps.

Fine Tuning

- Enters program mode starting with current program settings
1. Press reset button, confirmed by 1 beep.
 2. Within 2 seconds press reset button again, confirmed by 2 beeps.


Basic Reset

- Reverts to users programmed settings in normal operation mode.
1. Press reset button, confirmed by 1 beep.

Note... A basic reset should always be carried out after fitting batteries.

When would you carry out any of the above resets?


- No response from the limb when trying any of the knee activation sequences
- Limb mode unknown
- After changing or reconnecting batteries (Basic Reset)

Should unit require return to Blatchford, please remove and dispose of the batteries. (Disposal as per local waste handling regulations).

9 Alignment Advice and Stability Adjustment


The instructions in this section are for practitioner use only.

10 Fitting Advice



The instructions in this section are for practitioner use only.

10.1 Changing Stanceflex Rubbers





Static

Knee Flexion Following Heel Strike

10.2 Front Stop Adjustment - Adjustable Trunnion

The instructions in this section are for practitioner use only.

Front Stop Adjustment - Non-Adjustable Intelligent Prosthesis

(As fitted to older SmartIP models.)

No adjustment is necessary when using these cylinders. They require the front stop to be fully tightened as shown in Stage 1.

10.3 Further Adjustments

The instructions in this section are for practitioner use only.

Symptom	Solution
Knee does not stabilise	Increase sensitivity of ESK+
	Shift socket anteriorly
	Replace stabilizing mechanism if trailing edge of brake is damaged
	Check Stanceflex rubber and replace if necessary (See minimum thickness in compression in <i>Knee Flexion Following Heel Strike</i> section)
Knee does not release	Decrease sensitivity of ESK+
	Check correct adjustment of piston rod length
	Shift socket posteriorly
Heel rise too high during swing	See control cylinder IFU
Heel rise too low during swing	See control cylinder IFU
User experiences terminal impact	See control cylinder IFU
Device does not fully extend during swing	See control cylinder IFU



If after following these instructions you still have problems with the function please contact the sales team in your area for advice.

11 Assembly Instructions


The instructions in this section are for practitioner use only.

11.1 Changing Front Stop Buffer


When fitting a new front stop, always re-adjust, and if applicable check length of cylinder if fitted with an adjustable trunnion. (See *Front Stop Adjustment - Adjustable Trunnion.*)

11.2 Changing Shin Stops









1 Remove shin stops and remainder of adhesive. Clean and lightly abrade mating surfaces.

2 Apply Loctite 424 on mating surfaces. Press shin stops into place.









11.3 Removing ESK+ Mechanism

The instructions in this section are for practitioner use only.

<p>1 Remove all preload and record number of clicks</p> <p>Underside view of chassis</p> 	<p>2</p> 	<p>3</p> 
<p>4</p> 	<p>5</p> 	<p>6</p> 
<p>7</p>  <p>Oscillate during removal</p> <p>Brake band material may be wiped clean with paper, do not use solvents.</p>		<p>8</p>  <p>Clean bore with suitable solvent</p>

11.4 Replacement of ESK+ Mechanism

The instructions in this section are for practitioner use only.


<p>1</p>  <p>Offer up side plate to align slot in spindle (2) with dowel (1) but do not fit at this stage.</p>	<p>2</p>  <p>Screw spindle to chassis.</p>	
<p>3</p> <p>Oscillate during assembly</p>  <p>ORIENTATION IS IMPORTANT! Assemble flat sides and lettering (RH, LH) facing inwards.</p>	<p>4</p> 	<p>5</p> 
<p>6</p> 	<p>7</p> 	<p>8</p> <p>Re-apply preload with recorded number of clicks.</p> <p>Underside view of chassis.</p> 

12 Technical Data


Operating and Storage Temperature Range:	-10°C to 50°C (14°F to 122°F)
Component Weight:	
4-Bolt	1.15 kg (2 lb 8.6 oz)
Single Bolt Stanceflex	1.27 kg (2 lb 13 oz)
Activity Level:	2-3 (4)
Maximum User Weight:	
Activity Levels 2-3	125 kg (275 lb)
Activity Level 4	100 kg (220 lb 7.4 ozs)
Proximal Attachment:	4-Bolt or Single bolt Stanceflex (Blatchford)
Distal Alignment Attachment:	30 mm Tube Clamp or Male Pyramid (Blatchford)
Build Height:	(See diagrams below.)

Build Height

4-Bolt


Distal Pyramid


30 mm Distal Tube Clamp

Stanceflex

Distal Pyramid

30 mm Distal Tube Clamp

13 Ordering Information

	30 mm Distal Tube Clamp			Distal Pyramid		
4-bolt Adapter	SMARTIP	4B	30	SMARTIP	4B	PYR
	SmartIP Knee	4-bolt Adapter	Distal Tube Clamp 30 mm	SmartIP Knee	4-bolt Adapter	Distal Pyramid
Stanceflex	SMARTIP	SF	30	SMARTIP	SF	PYR
	SmartIP Knee	Stanceflex	Distal Tube Clamp 30 mm	SmartIP Knee	Stanceflex	Distal Pyramid

Replacement Parts

Item	Part. No.
Stanceflex Knee Assembly	239149
ESK Mechanism Alignment Pin	234914
Front Stop Pad	234903
ESK+ Stabilizing Mechanism Kit	239949
Shin Stop	233810

Stanceflex Rubber

Item	Part. No.
25° Pink (Soft)	231608
35° White	231606
45° Orange (Fitted as Standard)	231607
70° Blue (Hard)	231609

FCC Compliance Statement

This device complies with part 15 of the FCC Rules.

Operation is subject to the following two conditions:

1. This device may not cause harmful interference.
2. This device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Caution: Exposure to Radio Frequency Radiation. This device must not be co-located or operating in conjunction with any other antenna or transmitter.

Responsible party: Blatchford Inc.

1031 Byers Road, Miamisburg, OH 45342
800-548-3534 (toll free) | 937-291-3636
customerservice@blatchfordus.com

RSS Compliance Statement

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions:

1. this device may not cause interference, and
2. this device must accept any interference, including interference that may cause undesired operation of this device.

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Caution: Exposure to Radio Frequency Radiation.

The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website Blatchford Inc. 1031 Byers Road, Miamisburg, OH 45342.

800-548-3534 (toll free) | 937-291-3636 customerservice@blatchfordus.com

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: www.blatchford.co.uk

Medical Device

Single Patient – multiple use

Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a practitioner.

Warranty

This device is warranted for 24 months, the brake band for 12 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

See the Blatchford website for the current full warranty statement.

Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

This symbol indicates that the product contains electrical/electronic components and/or batteries that should not be disposed of in general waste or be incinerated at the end of the product's life.

At the end of the product's life, all electrical/electronic components and/or batteries should be recycled or disposed of in accordance with the current regulations for the handling of WEEE (Waste Electrical and Electronic Equipment), or equivalent local regulations. The remainder of the product should also be recycled where possible in accordance with local waste recycling regulations.

To help prevent potential harm to the environment or to human health from uncontrolled waste disposal, Blatchford offers a take-back service. Please contact Customer Services for details.

Retaining the Packaging Label

The practitioner is advised to keep the packaging label as a record of the device supplied.

Trademark Acknowledgements

SmartIP and Blatchford are registered trademarks of Blatchford Products Limited.

Manufacturer's Registered Address

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Blatchford Europe GmbH
Endolite India Ltd.