



# KX07

The rugged knee designed for outstanding stability on all terrain

### READY FOR ACTION

A day at the lake, a walk in the park or a hike in the woods, this durable design is good to go everywhere.

### STANDS UP TO NATURE

Whether the path is wet, muddy or dusty, simply rinse this waterproof knee and keep going.

### RESPONDS NATURALLY

The unique combination of the polycentric knee and hydraulic swing-and-stance provides optimum stability with smooth walking.

### PRECISION MOVEMENT

Moving in tight spaces, kneeling, getting in and out of chairs and cars—it's easier with the large flexion angle (160).

### INHERENT STABILITY

The easy initiation of flexion increases toe clearance during swing, reducing the risk of trips or falls with far less effort than a single axis knee.

### HIGH IMPACT ACTIVITIES

It withstands rough terrain and wet, muddy or dusty environments, making it ideal for gardening, hiking, camping, and mountain biking.

### KX07 matches brilliantly with these products

#### BLADEXT AND ELITE BLADE

BladeXT and Elite Blade provide lightweight and dynamic solutions for more active users.

#### HYDRAULIC ANKLES

Blatchford Hydraulic Technology provides a smooth, safe, and natural walking experience, helping to contribute to your long-term health.



Designed for what moves you  
[blatchfordus.com](http://blatchfordus.com)

Blatchford



# KX07


## On top of your world



Blatchford

### Technical Specification

- Activity Level 3-4
- Max. user weight 330lbs (K3) and 319lbs (K4)
- Build Height 8 5/8"
- Weight 2.9lb (excluding cover)
- Knee Flexion 160°
- Warranty 36 months



### Ordering Information

Part Number KX07



WALK WAYS

TRAIL WAYS

WATER WAYS



An annual visual inspection is recommended. Check for visual defects that may affect proper function. Maintenance must be carried out by competent personnel. Before carrying out any new activities of daily living, please check with your clinician whether specific training is required.

Blatchford

@blatchfordUS | blatchfordus.com

800 548 3534 | info@blatchfordus.com | Blatchford Inc., 1031 Byers Road, Miamisburg, Ohio 45342, USA.

779711849 Iss1 10/20. Information correct at time of print.

# KX07

All terrain living