



Sandal Toe

Product Range

Blatchford 

Foot Range Now With Sandal Toe

Natural look

Flexible footwear choice

Additional sizes*

Available in 2 skin tones



*Additional width sizes are now available for sizes 25, 26 and 27.

Blatchford's new Sandal Toe Footshell offers a more natural look and provides flexibility of footwear choice with its enhanced sandal toe design. It is available in both a Light and Dark shade to help better match skin tones and offers additional width options helping you to find the perfect fit.

Linx

Multi-award winning, integrated limb system, designed to deliver an experience that mimics the incredible and complex structure of the human leg.



| Product Code | Size | Side | Width* | Spring Set | Sandal Toe |
|--------------|------|------|--------|------------|------------|
| LINX | 25 | L | N | 3 | S |

Example: foot size 25, left, narrow, spring rating 3, sandal toe.



Echelon^{VAC}

Hydraulic ankle technology paired with an elevated vacuum system for a secure and comfortable socket connection.



| Product Code | Size | Side | Width* | Spring Set | Sandal Toe |
|--------------|------|------|--------|------------|------------|
| EVAC | 25 | L | N | 3 | S |

Example: foot size 25, left, narrow, spring rating 3, sandal toe.



Elan

Microprocessor controlled hydraulic ankle for terrain response and speed control.



| Product Code | Size | Side | Width* | Spring Set | Sandal Toe |
|--------------|------|------|--------|------------|------------|
| ELAN | 25 | L | N | 3 | S |

Example: foot size 25, left, narrow, spring rating 3, sandal toe.



Echelon

Biomimetic hydraulic ankle control and split toe for stability on uneven or sloping terrain.



| Product Code | Size | Side | Width* | Spring Set | Sandal Toe |
|--------------|------|------|--------|------------|------------|
| EC | 25 | L | N | 3 | S |

Example: foot size 25, left, narrow, spring rating 3, sandal toe.



Echelon^{VT}

Biomimetic hydraulic ankle with shock absorber, rotation and split toe for terrain conformance.

| Product Code | Size | Side | Width* | Spring Set | Axial Spring | Sandal Toe |
|--------------|------|------|--------|------------|--------------|------------|
| ECVT | 25 | L | N | 5 | 3 | S |

Example: foot size 25, left, narrow, spring rating 5, axial spring 3, sandal toe.



Javelin

Light weight blade, ideal for moderate impact sports and activities

| Product Code | Size | Side | Width* | Spring Set | Sandal Toe | + Top Housing |
|--------------|------|------|--------|------------|------------|---------------|
| JAV | 25 | L | N | 3 | S | Adapter |

Example: foot size 25, left, narrow, spring rating 3, sandal toe, adapter.



Esprit

Low profile and lightweight foot ideal for longer residual limbs.

| Product Code | Size | Side | Width* | Spring Set | Sandal Toe |
|--------------|------|------|--------|------------|------------|
| ESP | 25 | L | N | 3 | S |

Example: foot size 25, left, narrow, spring rating 3, sandal toe.



Epirus

Multi-axial foot with ankle joint for ground compliance and energy efficient response.

| Product Code | Size | Side | Width* | Spring Set | Sandal Toe |
|--------------|------|------|--------|------------|------------|
| EP | 25 | L | N | 3 | S |

Example: foot size 25, left, narrow, spring rating 3, sandal toe.



*Narrow (N) and Wide (W) available for sizes 25-27 only.
For dark tone add suffix D.

Footshell Width Guide

Use the below table to help choose the correct size when ordering, for sizes 25, 26 and 27 you can now choose between Narrow (N) and Wide (W) fittings.



| Size | 22 | 23 | 24 | 25 | 26 | 27 | 28 | 29 | 30 |
|-----------------|------|------|------|----------------------|----------------------|----------------------|-------|-------|-------|
| Footshell Width | 76mm | 78mm | 80mm | 82mm (N) 88mm (W) | 83mm (N) 91mm (W) | 91mm (N) 99mm (W) | 100mm | 102mm | 105mm |
| Old Footshell | 77mm | 77mm | 80mm | 82mm | 82mm | 87mm | 90mm | 95mm | 101mm |

NB: Widths measured at widest point.

Spring Kits and Footshell Information

If you wish to order spring kits or footshells separately, please refer to the order codes below to ensure you receive the correct item.

Spring Kits (for Linx, Elan, Echelon, EchelonVT, EchelonVAC, Esprit, Epirus):

| Foot Size | Foot Size | Foot Size | Foot Size |
|-----------------|-----------------|-----------------|-----------------|
| 22-24 | 25-26 | 27-28 | 29-30 |
| Set 1 - 539801S | Set 1 - 539810S | Set 1 - 539819S | Set 1 - 539828S |
| Set 2 - 539802S | Set 2 - 539811S | Set 2 - 539820S | Set 2 - 539829S |
| Set 3 - 539803S | Set 3 - 539812S | Set 3 - 539821S | Set 3 - 539830S |
| Set 4 - 539804S | Set 4 - 539813S | Set 4 - 539822S | Set 4 - 539831S |
| Set 5 - 539805S | Set 5 - 539814S | Set 5 - 539823S | Set 5 - 539832S |
| Set 6 - 539806S | Set 6 - 539815S | Set 6 - 539824S | Set 6 - 539833S |
| Set 7 - 539807S | Set 7 - 539816S | Set 7 - 539825S | Set 7 - 539834S |
| Set 8 - 539808S | Set 8 - 539817S | Set 8 - 539826S | Set 8 - 539835S |

(for Javelin Spring Kits):

| Foot Size | Foot Size | Foot Size | Foot Size |
|-----------------|-----------------|-----------------|-----------------|
| 22-24 | 25-26 | 27-28 | 29-30 |
| Set 1 - 539301S | Set 1 - 539510S | Set 1 - 539519S | Special |
| Set 2 - 539302S | Set 2 - 539511S | Set 2 - 539520S | Order Only |
| Set 3 - 539303S | Set 3 - 539512S | Set 3 - 539521S | Set 3 - 539530S |
| Set 4 - 539304S | Set 4 - 539513S | Set 4 - 539522S | Set 4 - 539531S |
| Set 5 - 539305S | Set 5 - 539514S | Set 5 - 539523S | Set 5 - 539532S |
| Set 6 - 539306S | Set 6 - 539515S | Set 6 - 539524S | Set 6 - 539533S |
| Set 7 - 539307S | Set 7 - 539516S | Set 7 - 539525S | Set 7 - 539534S |
| Set 8 - 539308S | Set 8 - 539517S | Set 8 - 539526S | Set 8 - 539535S |
| Set 9 - - | Set 9 - 539518S | Set 9 - 539527S | Set 9 - 539536S |

Footshells*:

| Foot size Left/Right | Narrow | Wide |
|----------------------|----------|----------|
| 22L | 539038S | - |
| 22R | 539039S | - |
| 23L | 539040S | - |
| 23R | 539041S | - |
| 24L | 539042S | - |
| 24R | 539043S | - |
| 25L | 539044SN | 539044SW |
| 25R | 539045SN | 539045SW |
| 26L | 539046SN | 539046SW |
| 26R | 539047SN | 539047SW |
| 27L | 539048SN | 539048SW |
| 27R | 539049SN | 539049SW |
| 28L | - | 539050S |
| 28R | - | 539051S |
| 29L | - | 539052S |
| 29R | - | 539053S |
| 30L | - | 539054S |
| 30R | - | 539055S |

*For dark tone add suffix D

+44 (0) 1256 31660 | customer.service@blatchford.co.uk
Blatchford Products Limited, Unit D Antura, Kingsland Business Park,
Basingstoke, RG24 8PZ, United Kingdom.

@blatchfordgrp | blatchford.co.uk

203266220 Iss1 09/19. Information correct at time of print.

Blatchford